[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |
[Руководство по эксплуатации](#) | [사용지침](#) | [Bruksanvisning](#) | [使用説明書](#) |
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

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Table of Contents

1 Introduction	43	16 Calibration service	78
1.1 Scope of delivery.....	43	17 Warranty	79
1.2 Terms of use.....	43	18 Disposal	79
2 Safety Instructions	44		
2.1 General safety instructions.....	44		
2.2 Function.....	44		
2.3 Limitations of Use	46		
2.4 Operating Limitations.....	46		
2.5 Operating Exclusions	46		
2.6 Storage Conditions	47		
2.7 Recommended Application Range..	47		
2.8 Dispenser Selection Chart	48		
3 Functional and operating elements	50		
4 Assembly	52		
4.1 First Steps.....	52		
4.2 Priming	53		
5 Operation	55		
5.1 Dispensing.....	55		
5.2 Accessories	56		
6 Error limits	58		
7 Checking the Volume (Calibration)	59		
8 Adjustment	60		
8.1 Digital model	61		
8.2 Analog-adjustable model	62		
9 Cleaning	62		
9.1 Cleaning, Analog-adjustable model, Fixed-volume model	63		
9.2 Cleaning Digital model.....	64		
9.3 Replacing discharge tube/valves	66		
10 Autoclaving	67		
11 Troubleshooting	68		
12 Product markings	69		
13 Ordering Information	70		
14 Accessories/spare parts	72		
15 Repairs	77		
15.1 Sending for repair	77		

1 Introduction

1.1 Scope of delivery

Dispensette® S or Dispensette® S Organic bottle-top dispenser, for GL 45 threaded bottles, dosing cannula or dosing cannula with back-dosing valve, telescopic aspiration tube, back-dosing tube (optional for devices with back-dosing valve), assembly wrench, various bottle adapters, a quality certificate and these instructions for use.

Nominal volume ml	Adapter for bottle thread, PP	Filling tube length mm
1, 2, 5, 10	GL 24-25, GL 28, S 28, GL 32-33, GL 38, S 40	125- 240
25 / 50, 100	GL 32-33, GL 38, S 40	170- 330

1.2 Terms of use

- Carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-do-date versions of the operating manual on our website: www.brand.de.

1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

1.2.2 Symbols

Symbol	Meaning
	Hazardous area

1.2.3 Format

Format	Meaning	Format	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

2 Safety Instructions

2.1 General safety instructions

Please read carefully!

The instrument Dispensette® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and observe this operating manual before using the device.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
3. Observe all specifications provided by reagent manufacturers.
4. When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.
5. Use the instrument only for dispensing liquids, with strict regard to the defined limitations of use and operating limitations. Comply with the operating exclusions (Operating Exclusions, p. 46)! If in doubt, contact the manufacturer or supplier.
6. Always perform work in a manner that does not endanger yourself or other people. When dispensing, the discharge tube must always point away from you or any other person. Avoid splattering. Only use suitable vessels.
7. Never press down the piston when the discharge tube closure is attached.
8. Never remove the discharge tube while the dispensing cylinder is filled.
9. Reagents can accumulate in the screw cap of the discharge tube. Thus, the screw cap should be cleaned regularly.
10. For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over.
11. Never carry the mounted instrument by the cylinder sleeve or the valve block. Breakage or loosening of the cylinder can lead to personal injury from chemicals (First Steps, p. 52, Fig. 3).
12. Never use force. Use smooth gentle movements to operate the piston upwards and downwards.
13. Use only original accessories and original replacement parts. Do not make any technical modifications. Do not dismantle the instrument any further than is described in the operating manual!
14. Always check that the instrument is in proper working condition before use. The user can come into contact with media if the instrument has been insufficiently cleaned or inspected. If there is a sign of a potential malfunction (e.g., piston difficult to move, sticking valves or leakage), immediately stop dispensing and consult the Troubleshooting section of this manual (Troubleshooting, p. 68). Contact the manufacturer, if necessary.

2.2 Function

With the Dispensette® S and Dispensette® S Organic bottle-top dispensers, liquids can be dispensed directly from the supply bottle. Available in Digital, Analog and Fixed models. The instruments are marked DE-M and optionally equipped with recirculation valve.

2.2.1 Dispensette S (red color code)

Digital

Analog

Fixed

2.2.2 Dispensette S Organic (yellow color code)

Digital

Analog

Fixed

2.2.3 Operation

NOTICE

Dispensing hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials:

Dispensette® S

Borosilicate glass, Al₂O₃-ceramic, ETFE, FEP, PFA, PTFE, platinum-iridium, PP (screw cap).

Dispensette® S Organic

Borosilicate glass, Al₂O₃ceramic, ETFE, FEP, PFA, PTFE, tantalum, PP (screw cap).

Alternatively ETFE/PTFE bottle adapters can be used (Accessories/spare parts, p. 72). The suitability of ETFE/PTFE bottle adapters must be checked by the user.

2.3 Limitations of Use

This instrument is designed for dispensing liquids, observing the following physical limits:

- Operating temperature from +15 °C to +40 °C (from 59 °F to 104 °F) of instrument and reagent
- Vapor pressure up to max. 600 mbar. Aspirate slowly above 300 mbar, in order to prevent the liquid from boiling
- Kinematic viscosity up to 500 mm²/s (dynamic viscosity [mPas] = kinematic viscosity [mm²/s] x density [g/cm³])
- Density up to 2.2 g/cm³

2.4 Operating Limitations

- Liquids that form deposits (e.g., crystallizing solutions or concentrated alkaline solutions) may make the piston difficult to move or may cause jamming. If the piston movement becomes sluggish or stiff, the instrument should be cleaned immediately. See also Cleaning, p. 62.
- When dispensing inflammable media, make sure to avoid to buildup of static charge (e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth).
- If hazardous media (e.g., hydrofluoric acid or bromine) are dispensed, replace the dispensing unit after approx. 3000 complete strokes. Replacement may also be necessary earlier depending on the medium dispensed and the frequency of use. See also Replacing the dispensing cartridge.
- The instrument is designed for general laboratory applications and complies with the requirements of the relevant standards (e.g., DIN EN ISO 8655). Compatibility of the instrument for a specific application (e.g., trace material analysis or food sector) must be checked by the user. Approvals for specific applications (e.g., for the production and administration of food, pharmaceuticals, or cosmetics) have not been granted.

2.5 Operating Exclusions

2.5.1 Dispensette S

Never use Dispensette® S for:

- Liquids that attack Al₂O₃-ceramic, ETFE, FEP, PFA, and PTFE (e.g., dissolved sodium azide*)
- Liquids that attack borosilicate glass (e.g., hydrofluoric acid)
- Liquids that are decomposed catalytically by platinum-iridium (e.g., H₂O₂)
- Hydrochloric acid > 20% and nitric acid > 30%
- Tetrahydrofuran
- Trifluoroacetic acid
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., activated charcoal) because solid particles may clog or damage the instrument
- Liquids that attack PP (screw cap and bottle adapter)**.

* Dissolved sodium azide permitted up to a concentration of max. 0.1 %.

** Alternatively, ETFE/PTFE bottle adapters and ETFE/PTFE screw caps can be used (accessories from Accessories/spare parts, p. 72). The suitability of ETFE/PTFE bottle adapters and ETFE/PTFE screw caps must be checked by the user.

2.5.2 Dispensette S Organic

Never use Dispensette® S Organic for:

- Liquids that attack Al_2O_3 -ceramic, tantalum, ETFE, FEP, PFA, and PTFE (e.g., dissolved sodium azide*)
- Liquids that attack borosilicate glass (e.g., hydrofluoric acid)
- Bases and salt solutions
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., activated charcoal) because solid particles may clog or damage the instrument
- Liquids that attack PP (screw cap and bottle adapter)**.

* Dissolved sodium azide permitted up to a concentration of max. 0.1 %.

** Alternatively, ETFE/PTFE bottle adapters and ETFE/PTFE screw caps can be used (accessories from Accessories/spare parts, p. 72). The suitability of ETFE/PTFE bottle adapters and ETFE/PTFE screw caps must be checked by the user.

2.6 Storage Conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from -20 to $50^\circ C$ (-4 to $122^\circ F$).

2.7 Recommended Application Range

Dispensette® S offers a wide range of applications for dispensing aggressive reagents (e.g., concentrated acids such as H_3PO_4 and H_2SO_4 , bases such as NaOH and KOH, and salt solutions as well as various polar solvents).

Dispensette® S Organic is used for dispensing organic solvents (e.g., chlorinated and fluorinated hydrocarbons such as trichlorotrifluoroethane and dichloromethane or acids such as concentrated HCl and HNO_3 as well as for trifluoroacetic acid (TFA), tetrahydrofuran (THF) and peroxides).

NOTICE

Dispenser selection

For guidelines on selecting the right dispenser, please observe the corresponding exclusions of use and the following 'Dispenser Selection Chart'.

NOTICE

Dispensing hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

2.8 Dispenser Selection Chart

NOTICE

Hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the Dispensette S Trace Analysis bottle-top dispenser with platinum-iridium valve spring. See separate instructions for use at www.brand.de/om.

	Dispensette® S	Dispensette® S Organic
Acetaldehyde	✓	✓
Acetone	✓	✓
Acetonitrile	✓	✓
Acetophenone	—	✓
Acetylacetone	✓	✓
Acetyl chloride	—	✓
Acrylonitrile	✓	✓
Acrylic acid	✓	✓
Adipic acid	✓	—
Allyl alcohol	✓	✓
Aluminum chloride	✓	—
Formic acid, ≤ 100%	—	✓
Amino acids	✓	—
Ammonia, ≤ 20%	✓	✓
Ammonia, 20–30%	—	✓
Ammonium chloride	✓	—
Ammonium fluoride	✓	—
Ammonium sulfate	✓	—
n-amyl acetate	✓	✓
Amyl alcohol (pentanol)	✓	✓
Amyl chloride (chloropentane)	—	✓
Aniline	✓	✓
Essential oils (aromatic oils)	—	✓
Barium chloride	✓	—
Benzaldehyde	✓	✓
Gasoline (petroleum spirit), boiling point 70–180°C	—	✓
Benzoic acid methyl ester	✓	✓
Benzene	✓	✓
Benzoyl chloride	✓	✓
Benzyl alcohol	✓	✓
Benzylamine	✓	✓
Benzyl chloride	✓	✓
Boric acid, ≤ 10%	✓	✓
Pyruvic acid	✓	✓
Bromobenzene	✓	✓
Bromonaphthalene	✓	✓
Hydrobromic acid	—	✓
Butanediol	✓	✓
1-butanol	✓	✓
Butyric acid	✓	✓
Butyl acetate	✓	✓

	Dispensette® S	Dispensette® S Organic
Butylamine	✓	✓
Butyl methyl ether	✓	✓
Calcium carbonate	✓	—
Calcium chloride	✓	—
Calcium hydroxide	✓	—
Calcium hypochlorite	✓	—
Chloroacetaldehyde, ≤ 45%	✓	✓
Chloroacetone	✓	✓
Chlorobenzene	✓	✓
Chlorobutane	✓	✓
Chloroacetic acid	✓	✓
Chloronaphthalene	✓	✓
Chloroform	—	✓
Chlorosulfonic acid	—	✓
Chromic acid, ≤ 50%	✓	✓
Chromosulfuric acid	✓	—
Cumene (isopropylbenzene)	✓	✓
Cyclohexane	—	✓
Cyclohexanone	✓	✓
Cyclopentane	—	✓
Decane	✓	✓
1-decanol	✓	✓
Dibenzyl ether	✓	✓
Dichlorobenzene	✓	✓
Dichloroacetic acid	—	✓
Dichloroethane	—	✓
Dichloroethylene	—	✓
Dichloromethane	—	✓
Diesel oil (fuel oil), boiling point 250–350°C	—	✓
Diethanolamine	✓	✓
Diethylamine	✓	✓
1,2-diethylbenzene	✓	✓
Diethylene glycol	✓	✓
Diethyl ether	—	✓
Dimethylaniline	✓	—
Dimethylformamide (DMF)	✓	✓
Dimethyl sulfoxide (DMSO)	✓	✓
1,4-dioxane	—	✓
Diphenyl ether	✓	✓
Acetic acid, ≤ 96%	✓	✓
Acetic acid, 100% (glacial acetic acid)	✓	✓
Acetic anhydride	—	✓
Ethanol	✓	✓
Ethanolamine	✓	✓
Ethyl acetate	✓	✓
Ethylbenzene	—	✓
Ethylene chloride	—	✓
Ethyl methyl ketone	✓	✓
Fluoroacetic acid	—	✓
Formaldehyde, ≤ 40%	✓	—
Formamide	✓	✓
Glycol (ethylene glycol)	✓	✓
Glycolic acid, ≤ 50%	✓	—
Glycerin	✓	✓
Urea	✓	—

	Dispensette® S	Dispensette® S Organic
Fuel oil (diesel oil), boiling point 250–350°C	—	✓
Heptane	—	✓
Hexane	—	✓
Hexanol	✓	✓
Hexanoic acid	✓	✓
Hydroiodic acid, ≤ 57% **	✓	✓
Isoamyl alcohol	✓	✓
Isobutanol	✓	✓
Isooctane	—	✓
Isopropanol (2-propanol)	✓	✓
Isopropyl ether	✓	✓
Potassium chloride	✓	—
Potassium dichromate	✓	—
Potassium hydroxide	✓	—
Potassium permanganate	✓	—
Cresol	—	✓
Copper sulfate	✓	—
Methanol	✓	✓
Methoxybenzene	✓	✓
Methyl butyl ether	✓	✓
Methylene chloride	—	✓
Methyl formate	✓	✓
Methyl propyl ketone	✓	✓
Lactic acid	✓	—
Mineral oil (engine oil)	✓	✓
Monochloroacetic acid	✓	✓
Sodium acetate	✓	—
Sodium chloride	✓	—
Sodium dichromate	✓	—
Sodium fluoride	✓	—
Sodium hypochlorite	✓	—
Sodium hydroxide, ≤ 30%	✓	—
Nitrobenzene	✓	✓
Oleic acid	✓	✓
Oxalic acid	✓	—
Pentane	—	✓
Perchloroethylene	—	✓
Perchloric acid	✓	✓
Peracetic acid	—	✓
Petroleum ether, boiling point 40–70°C	—	✓
Petroleum, boiling point 180–220°C	—	✓
Phenol	✓	✓
Phenylethanol	✓	✓
Phenylhydrazine	✓	✓
Phosphoric acid, 85	✓	✓
Phosphoric acid, 85	✓	✓
Piperidine	✓	✓
Propionic acid	✓	✓
Propylene glycol (propanediol)	✓	✓
Pyridine	✓	✓

	Dispensette® S	Dispensette® S Organic
Salicylaldehyde	✓	✓
Nitric acid, ≤ 30%	✓	✓
Nitric acid, 30–70% */**	—	✓
Hydrochloric acid, ≤ 20%	✓	✓
Hydrochloric acid, 20–37% **	—	✓
Sulfuric acid, ≤ 98%	✓	✓
Sulfuric acid, 98%, 1:1	✓	✓
Silver acetate	✓	—
Silver nitrate	✓	—
Scintillation cocktail	✓	✓
Turpentine	—	✓
Tetrachloroethylene	—	✓
Carbon tetrachloride	—	✓
Tetrahydrofuran (THF) */**	—	✓
Tetramethylammonium hydroxide	✓	—
Toluene	—	✓
Trichlorobenzene	—	✓
Trichloroacetic acid	—	✓
Trichloroethane	—	✓
Trichloroethylene	—	✓
Trichlorotrifluoroethane	—	✓
Triethanolamine	✓	✓
Triethylene glycol	✓	✓
Trifluoroacetic acid (TFA)	—	✓
Trifluoroethane	—	✓
Hydrogen peroxide, ≤ 35%	—	✓
Tartaric acid	✓	—
Xylene	—	✓
Zinc chloride, ≤ 10%	✓	—
Zinc sulfate, ≤ 10%	✓	—

Legend:

✓ = The device is suitable for the medium

— = The device is not suitable for the medium

This table has been carefully tested and is based on the most current information available. Always observe the operating manual of the instrument and the specifications provided by the reagent manufacturer. In addition to the chemicals listed above, a number of organic or inorganic saline solutions (e.g. biological buffers), biological detergents as well as media for cell cultures can be dispensed. If you need chemical hazard statements that are not on the list, you are welcome to contact BRAND.

Last updated: 0522/14

* Use ETFE/PTFE bottle adapter

** Use PTFE sealing ring for valve block


3 Functional and operating elements

Dispensette® S Digital

- 1 Volume adjustment
- 2 Cover plate
- 3 Housing
- 4 Protective sleeve/dispensing cylinder
- 5 Piston unit
- 6 Valve block
- 7 Valve block adapter (GL 45 bottle thread)
- 8 Recirculation tube
- 9 Telescopic filling tube
- 10 Latch
- 11 Safety lock
- 12 Piston mounting nut
- 13 Knob
- 14 Screw cap
- 15 Discharge tube
- 16 Valve lever, recirculation valve

Dispensette® S Analog

- 1 Cover plate
- 2 Piston seat
- 3 Housing
- 4 Pointer
- 5 Volume adjustment
- 6 Protective sleeve/dispensing cylinder
- 7 Valve block
- 8 Valve block adapter (GL 45 bottle thread)
- 9 Recirculation tube
- 10 Telescopic filling tube
- 11 Screw cap
- 12 Discharge tube
- 13 Valve lever, recirculation valve
- 14 Piston unit

Telescoping filling tube and recirculation tube

Assembly tool

4 Assembly

4.1 First Steps

⚠ WARNING

Follow the safety instructions

- > Wear protective clothing, eye protection and protective gloves!
- > Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- > Follow all safety instructions and comply with the limitations of use, see Limitations of Use, p. 46.
- > Comply with the operating limitations, see Operating limitations, p. 46.

1. Mounting the filling tube/recirculation tube

- a. Adjust the length of the telescoping filling tube to the bottle height and attach it.
 - b. Center and attach the filling tube (part with smaller diameter) carefully to avoid damaging the nozzle.
- ⇨ If a discharge tube with a recirculation valve is used, the recirculation tube must also be installed.
- c. Insert the recirculation tube with the opening pointing outward.

2. Mount the instrument on a bottle and align

NOTICE

To avoid tipping over, use a bottle stand for small bottles.

- a. Screw the instrument (GL 45 threads) onto the reagent bottle, and then align the discharge tube with the bottle label. This is done by rotating the valve block with the discharge tube.

3. Transporting the instrument

NOTICE

Polypropylene (PP) bottle adapters are supplied with the instrument. These may be used only for media that do not attack the PP. Alternatively ETFE/PTFE bottle adapters can be used (Accessories/spare parts, p. 72). The suitability of ETFE/PTFE bottle adapters must be checked by the user.

- a. For bottles with other thread sizes, select a suitable bottle adapter.
- b. Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- c. When mounted to a reagent bottle, always carry the instrument as shown in the figure!

4.2 Priming

⚠ WARNING

Follow during every use, especially with hazardous media

- Wear protective clothing, eye protection and appropriate hand protection!
- Never press down the piston when the screw cap is screwed on!
- Avoid splashing the reagent!
- To avoid splashes dispense slowly.
- Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use, p. 46 and Operating exclusions, p. 46.

NOTICE

Before using the instrument for the first time, ensure it is rinsed carefully and discard the first few samples dispensed. Avoid splashes, as ethanol and glycerin residues may be present in the instrument.

Instruments with recirculation valve

- a. Open the screw cap of the dispensing tube.

- b. Turn valve to 'Recirculate'.

- c. For deaerating gently pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this process at least 5 times.

- d. Turn valve to 'Dispense'.

- e. To avoid splashes when deaerating hold the discharge tube on the inner wall of a suitable receiving vessel and dispense liquid to vent deaerate the discharge tube until it is bubble-free. Wipe away any remaining drops from the discharge tube.

Instruments without recirculation valve

- a. Open the screw cap of the discharge tube (see 'instrument with recirculation valve', Fig. a.). To avoid splashes, hold discharge tube orifice on the inner wall of a suitable receiving vessel.
- b. For deaerating pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this procedure approximately 5 times until the discharge tube is bubble-free.

5 Operation

5.1 Dispensing

1. Selecting the volume

Digital: Rotate the volume-setting wheel until the desired volume is indicated (mechanical counter).

Analog: Loosen the volume selector thumb screw $\frac{3}{4}$ turn (1), set the pointer to the desired volume (2) and then retighten the volume thumb screw (3).

Fixed: The volume is non-adjustable and cannot be changed.

2. Dispensing

⚠ WARNING

Follow during every use, especially with hazardous media

- > Wear protective clothing, eye protection and appropriate hand protection!
- > Never press down the piston when the screw cap is screwed on!
- > Avoid splashing the reagent!
- > To avoid splashes dispense slowly.
- > Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- > Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use, p. 46 and Operating exclusions, p. 46.

- a. Unscrew the screw cap of the dispensing tube.
- b. When using instruments equipped with a recirculation valve, turn the valve to 'Dispensing'.
- c. Hold the discharge tube orifice on the inner wall of a suitable receiving vessel.

- d. Gently lift the piston until the upper stop and then depress piston slowly and steadily with minimal force until the lower stop.
- e. Wipe off the discharge tube against the inner wall of the receiving vessel.
- f. Seal the dispensing tube with the screw cap.

NOTICE

After using the piston, always press it down to the lower stop (parking position). If the piston is not pressed down to the lower stop, unintentional media leaks may occur.

NOTICE

The filled status of the instrument must be specially marked during cleaning.

5.2 Accessories

5.2.1 Flexible discharge tube with recirculation valve

For serial dispensing, the flexible discharge tubing for the bottle-top dispensers Dispensette® S and Dispensette® S Organic can also be used (Accessories/spare parts, p. 72).

The specified accuracy and coefficient of variation of the instrument are only obtained for volumes > 2 ml and by gently approaching the upper and lower stops. The coil of the tubing can be stretched to a maximum length of 800 mm. Before use, the tube must lie in regular loops and must not be twisted. The operating exclusions for the respective instrument in use apply.

Assembly



⚠ WARNING

Use only undamaged tubing

There should be no visible damage to the discharge tube (e.g. kinks, etc.). Each time you are going to use the tubing, examine it carefully.

- > To dispense aggressive liquids, you should take safety measures in addition to the normal precautions. We recommend use of a protective shield.
- > The bottle must be supported using a bottle stand.
- > To help avoid reagent splashing from the tube, always grip the tube firmly by the handle and replace into the holder after use.
- > For cleaning rinse the tube carefully.
- > Do not dismantle!

Prerequisite:

- If the instrument was in use, it must be cleaned before installing the flexible discharge tube (Cleaning, p. 62).
- a. Instruments with a recirculating valve should be set to 'Recirculate', and the valve lever pulled upwards to remove it.
- b. Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.
- c. Push the flexible discharge tube holder from the bottom of the valve block and tighten it. In doing so, the instrument must not be mounted to a bottle. Install the receiver tube.
- d. Press the plug of the recirculation valve downwards.
- e. Slide the flexible discharge tube housing into the valve block up to the stop.
- f. Push the housing all the way down.
- g. Put on the valve toggle made for the discharge valve and press it in firmly. In doing so, pay attention to the color coding and label (see assembly instructions 'Flexible discharge tubing for Dispensette® S').

NOTICE

Use a bottle stand (Accessories/spare parts, p. 72).

5.2.2 Drying tube

Use of a drying tube, filled with a suitable absorbent (purchased separately), might be necessary for moisture- and CO₂- sensitive media.

(Accessories/spare parts, p. 72)

Assembly

- a. Use the assembly tool to unscrew the air vent cap.
- b. Screw in the filled drying tube.

c.

- c. Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

NOTICE

If necessary, seal the threads of the drying tube, the bottle and/or the bottle adapter with PTFE tape.

5.2.3 Sealing ring for valve block

For highly volatile media we recommend to seal the connection from valve block to bottle with the PTFE sealing ring and PTFE tape (Accessories/spare parts, p. 72).

Assembly

a.

- a. Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

5.2.4 Air vent cap for micro filter with Luer-cone

For sterile media we recommend the air vent cap with Luer-cone to attach a micro filter. This provides increased protection against contamination by displacement air (Accessories/spare parts, p. 72).

Assembly

a.


- a. Unscrew the air vent cap (see 'Assembly Drying tube', Fig. a.).
- b. Screw in the air vent cap with a Luer cone.
- c. Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.
- d. Insert a commercially available sterile filter into the Luer cone.

c.

6 Error limits

Error limits in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20 °C/68 °F) of instrument, ambient environment and distilled water. Testing takes place according to DIN EN ISO 8655-

6 with a completely filled instrument and with uniform and smooth dispensing.

Error limits

Nominal volume ml	A* $\leq \pm$ %	μl	CV* \leq %	μl
1	0.5	5	0.1	1
2	0.5	10	0.1	2
5	0.5	25	0.1	5
10	0.5	50	0.1	10
25	0.5	125	0.1	25
50	0.5	250	0.1	50
100	0.5	500	0.1	100

$$A_T = \frac{V_N}{V_T} \cdot A_N$$

Partial volume

The percentage values for A and CV are relative to the nominal volume (V_N) and must be converted for partial volumes (V_T).

i.e.	Volume	A* $\leq \pm$ %	μl	CV* \leq %	μl
V_N	25.0	0.5	125	0.1	25
$V_p = 50\% N$	12.5	1.0	125	0.2	25
$V_p = 10\% N$	2.5	5.0	125	1.0	25

*A = Accuracy, CV = Coefficient of Variation

NOTICE

The error limits in DIN EN ISO 8655-5 are significantly lower. The maximum error for a single measurement can be approximated from the sum of error limits $EL = A + 2 \times CV$ (for 25 ml: $125 \mu\text{l} + 2 \times 125 \mu\text{l} = 175 \mu\text{l}$).

7 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. This time frame should be adjusted to correspond with individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de. For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. Demo version can be downloaded from www.brand.de. Gravimetric volume testing according to DIN EN ISO 8655-6 (for measurement conditions, see Error limits, p. 58) is performed as follows:

1. Prepare the instrument

Clean the instrument (Cleaning, p. 62), fill it with distilled H_2O and then prime it carefully.

2. Check the volume

- At 10 pipetting series with distilled H₂O and weighings in 3 volume ranges (100 %, 50 %, 10 %) are recommended
- For discharge, depress the piston slowly and steadily without force until the lower stop
- Wipe off the tip of titration tube.
- Weigh the pipetted amount with an analysis scale. (Please refer to the user manual of the scale manufacturer.)
- Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Coefficient of variation*:

$$CV\% = \frac{100 s}{\bar{V}}$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

NOTICE

Test instructions (SOPs) are available for download at www.brand.de.

8 Adjustment

After a long period of usage, an adjustment of the instrument might be necessary.

- Calibrate for example at nominal volume (Checking the Volume (Calibration)).
- Calculate mean volume (actual value) (Checking the Volume (Calibration)).
- Adjust the instrument (to the calculated mean volume).
- After the adjustment, further calibration is necessary to confirm the appropriate adjustment.

Example

The gravimetric check gives an actual value of 9.90 ml for a 10 ml instrument set for a nominal volume of 10.00 ml.


8.1 Digital model

- a.** Slide the latch to the left and remove the front housing shell.

- b.** Pull out the safety lock. The adjustment cover will then come off. Discard the adjustment cover.

- c.** Pull the red knob to disengage the gears. Set the display to actual delivered volume (e.g., 9.90 ml).

- d.** Press in the red knob and then the safety lock again.

- e.** Replace housing and slide the latch to the right.
 ⇒ Adjustment is completed.
 ⇒ Alteration of the factory setting is indicated by a red recalibration flag (circle).

Adjustment range

Nominal volume [ml]	Digital max. +/- [μ l]
1	12
2	24
5	60
10	120
25	300
50	600
100	–

8.2 Analog-adjustable model

- a. Insert the pin of the assembly tool into the cover plate (Pos. 1), and break it off with a rotating motion. Discard the adjustment cover.
- b. Insert the pin of the assembly tool into the adjustment screw and rotate to the left in order to increase the dispensing volume, or rotate to the right to decrease the dispensing volume (e.g. for an actual value of 9.97 ml, rotate approx. 1/2 turn to the left).
- c. Adjustment is completed.
 - ⇒ The change in the adjustment is indicated by a red disk (circle in figure).

Adjustment range

Nominal volume	Analog-adjustable/Fixed-volume max. +/- [μ l]	One rotation corresponds to [μ l]
1 ml	6	~ 8
2 ml	12	~ 16
5 ml	30	~ 40
10 ml	60	~ 80
25 ml	150	~ 130
50 ml	300	~ 265
100 ml	600	~ 400

9 Cleaning

⚠ WARNING

Components filled with reagent

The cylinder, valves, telescopic filling tube and discharge tube contain reagent!

- > Never remove the discharge tube while the dispensing cylinder is filled.
- > Point the valves and tube openings away from your body.
- > Wear protective clothing, eye protection and protective gloves!

The instrument must be cleaned in the following situations to assure correct operation:

- immediately when the piston is difficult to move
- prior to a reagent change
- prior to long term storage
- prior to changing the valve
- regularly when using liquids which form deposits (e.g., crystallizing liquids)

- prior to dismantling the instrument
- regularly when liquids accumulate in the screw cap.
- prior to autoclaving

For proper cleaning and removal of any deposits in the parts through which liquids pass, also always completely with draw the piston from the cylinder after rinsing with a suitable cleaning solution. If necessary, the parts can also be cleaned in an ultrasonic bath.

9.1 Cleaning, Analog-adjustable model, Fixed-volume model

1. Completely emptying the instrument

- Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in both the 'dispense' and 'recirculate' settings.

2. Rinsing the instrument

- Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.

3. Removing pistons


NOTICE

Avoid swapping pistons

The pistons and cylinders are individually matched, and should not be swapped with pistons from other instruments!

- Hold the housing securely and unscrew the piston seat completely by turning it to the left.

- Pull out the piston carefully.

4. Cleaning the piston and cylinder and reassembling the instrument

NOTICE

Special instructions for Dispensette® S Organic

For the Dispensette® S Organic, slide the piston into the cylinder vertically with a rotating motion.

- Clean the piston and cylinder and reassemble the instrument. If necessary carefully remove deposits at the edge of the glass cylinder.
- Rinse the piston and cylinder with deionized water, and dry them carefully.
- Insert the piston completely into the cylinder and then reassemble the instrument. The piston can only be inserted into the cylinder when the discharge tube is attached.

9.2 Cleaning Digital model

1. Completely emptying the instrument

- Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in both the 'dispense' and 'recirculate' settings.

2. Rinsing the instrument

- Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.

3. Removing pistons


NOTICE

Avoid swapping pistons

The pistons and cylinders are individually matched, and should not be swapped with pistons from other instruments!

- Carry out assembly and disassembly at the maximum volume setting only. To do this, turn the volume adjustment wheel as far as it will go in + direction (= maximum volume).
- Slide the latch to the left.
- Remove the front housing shell.

- d. Use the assembly tool to loosen the piston mounting nut and pull the rear housing section backwards.
- e. Then carefully pull the piston out of the cylinder.

4. Cleaning the piston and cylinder and reassembling the instrument

NOTICE

Special instructions for Dispensette® S Organic



For the Dispensette® S Organic, slide the piston into the cylinder vertically with a rotating motion.

- a. Clean the piston and cylinder and reassemble the instrument. If necessary carefully remove deposits at the edge of the glass cylinder.
- b. Rinse the piston and cylinder with deionized water, and dry them carefully.
- c. Insert the piston completely into the cylinder and then reassemble the instrument. The piston can only be inserted into the cylinder when the discharge tube is attached.

NOTICE

Cleaning Digital model


The stop segment must engage underneath the cylinder's stop ring. When tightening the piston mounting nut using the assembly tool, use your thumbs to press the piston/cylinder unit back towards the back side of the housing.

9.3 Replacing discharge tube/valves

NOTICE

A functional test must always be carried out after exchanging components.

9.3.1 Replacing discharge tube

- Instruments with a recirculating valve should be set to 'Recirculate', and the valve lever pulled upwards to remove it.
- Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.
- Hold coupling piece of the new discharge tube and pull housing up. Push housing into the valve block until it meets the stop.
- Slide the discharge tube housing all the way down.
- For instruments with a recirculation valve, pull up the valve lever to the 'Recirculate' position, and press it in tightly.

9.3.2 Replacing valves

NOTICE

Always install the valve intended for the particular instrument model and size! (See ordering data, Ordering Information, p. 70) Identical filling valves but different discharge valves are used for Dispensette® S and Dispensette® S Organic. For easy identification, the discharge valves for the Dispensette® S Organic are marked with 'ORG'.

9.3.2.1 Replacing the discharge valve

NOTICE

The discharge valve has a safety ball that closes automatically when the discharge tube is removed. Attaching the discharge tube reopens the safety ball.

- After disassembling the discharge tube (see Replacing discharge tube, p. 66), use the assembly tool to unscrew the discharge valve.
- Firmly screw in the new discharge valve first by hand and then tighten it with the assembly tool. The thread should not be visible.

9.3.2.2 Replacing the filling valve

- a. Pull out the recirculation tube and the telescopic filling tube.
- b. Use the assembly tool to unscrew the filling valve.
- c. Screw in the new filling valve first by hand and then tighten it with the assembly tool.

9.3.2.3 Dislodging a stick valve ball

If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the valve ball is stuck.

In this case, dislodge the valve ball using light pressure; for example, with a 200 µl plastic pipette tip (see the figure at the side).

10 Autoclaving

The device can be autoclaved at 121°C (250°F), 2 bar, and a holding time of at least 15 minutes according to DIN EN 285.

- a. The instrument must be carefully cleaned prior to autoclaving (Cleaning, p. 62).
- b. Open the closure cap on the discharge tube, and for instruments with a recirculation valve, set the valve to 'Dispense'.
- c. Check that the filling valve is securely seated (Fig. 1a). In the Digital model, additionally check that the piston mounting is securely seated (Fig. 1b).
- d. To ensure unhindered access for the steam and to prevent the valve ball in the filling valve from possibly becoming stuck, hold the instrument with the discharge piston pressed vertically downward, and gently tap against the casing with your hand (Fig. 2). Then lay it horizontally in the autoclave. Be sure to avoid the instrument coming into contact with metal surfaces in the autoclave!

2

NOTICE

Do not reassemble the instrument until it has cooled down to room temperature (Cooling time approx. 2 hours). After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them. The effectiveness of autoclaving must be verified by the user.

11 Troubleshooting

Problem	Possible cause	Corrective action
Piston sluggish or stuck	Crystal deposits, Contaminants	Stop dispensing immediately. Loosen piston with circular motion, but do not disassemble. Carry out cleaning (Cleaning, p. 62).
Filling not possible	Volume adjusted to minimum setting	Set the desired volume (Dispensing, p. 55).
	Filling valve stuck	Unscrew the filling valve from the valve block and clean it. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it (Dislodging a stick valve ball, p. 67). If necessary, replace the filling valve.
Dispensing not possible	Discharge valve stuck	Unscrew the discharge valve from the valve block, clean it and replace if necessary. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it.
Discharge tube (with or without recirculation valve) cannot be mounted sufficiently	Discharge valve is not screwed in deeply enough	Tighten the discharge valve with the assembly tool until it meets the stop so that the threads are no longer visible.
Air bubbles in the instrument	Reagent with high vapor pressure has been drawn in too quickly	Slowly draw in reagent.
	Valve screw connections loose	Tighten the valves firmly with the assembly tool.
	The instrument has not been primed	Prime the instrument (Priming, p. 53).
	Filling tube is loose or damaged	Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and reconnect it or replace filling tube.
	Valves not firmly connected or damaged	Carry out cleaning (Cleaning, p. 62). Tighten the valves using the assembly tool.
Dispensed volume too low	Filling tube is loose or damaged	Carry out cleaning (Cleaning, p. 62). Push the filling tube on

Problem	Possible cause	Corrective action
		firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube (Replacing discharge tube, p. 66).
	Filling valve is loose or damaged	Carry out cleaning (Cleaning, p. 62). Tighten the valves using the assembly tool. If necessary, replace filling valves.
Leaking liquid between instrument and bottle	Recirculation tube not connected	Connect recirculation tube (First Steps, p. 52, Fig. 3).
	Volatile reagent dispensed without sealing ring	Install sealing ring for valve block (First Steps, p. 52).

12 Product markings

Symbol or number	Meaning
	General warning sign
	Note user manual
	Note eye protection
	Use hand protection
	Use protective clothing
XXZXXXXX	Serial number
	The device is marked in accordance with the German Weights and Measures Act and the Weights and Measures Ordinance. Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was added.
www.brand.de/ip	Patent information

13 Ordering Information

Dispensette® S, Digital

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.1 – 1	0.005	4600310	4600311
0.2 – 2	0.01	4600320	4600321
0.5 – 5	0.02	4600330	4600331
1 – 10	0.05	4600340	4600341
2.5 - 25	0.1	4600350	4600351
5 – 50	0.2	4600360	4600361

Dispensette® S, Analog-adjustable

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.1 - 1	0.02	4600100	4600101
0.2 - 2	0.05	4600120	4600121
0.5 - 5	0.1	4600130	4600131
1 - 10	0.2	4600140	4600141
2.5 - 25	0.5	4600150	4600151
5 - 50	1.0	4600160	4600161
10 - 100	1.0	4600170	4600171

Dispensette® S, Fixed-volume

Volume ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
1	4600210	4600211
2	4600220	4600221
5	4600230	4600231
10	4600240	4600 241
Fixed volumes by request: 0.5–100 ml (please specify when ordering!)	4600290	4600 291

NOTICE

Scope of supply, see Scope of delivery, p. 43

Dispensette® S Organic, Digital

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.5 - 5	0.02	4630330	4630331
1 - 10	0.05	4630340	4630341
2.5 - 25	0.1	4630350	4630351
5 - 50	0.2	4630360	4630361

Dispensette® S Organic, Analog-adjustable

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.5 - 5	0.1	4630130	4630131
1 - 10	0.2	4630140	4630141
2.5 - 25	0.5	4630150	4630151
5 - 50	1.0	4630160	4630161
10 - 100	1.0	4630170	4630171

Dispensette® S Organic, Fixed-volume

Volume ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
5	4630230	4630231
10	4630240	4630241
Fixed volumes by request: 2-100 ml (please specify when ordering!)	4630290	4630291

NOTICE**Dispensing hydrofluoric acid (HF)**

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

14 Accessories/spare parts

Bottle adapter

External thread	for bottle thread/ground joint	Material	Order No.
GL 32	GL 25	PP	704325
GL 32	GL 28/ S 28	PP	704328
GL 32	GL 30	PP	704330
GL 32	GL 45	PP	704345
GL 45	GL 32	PP	704396
GL 45	GL 35	PP	704431
GL 45	GL 38	PP	704397
GL 45	S* 40	PP	704343
GL 45	S* 50	PP	704350
GL 45	S* 54	PP	704430
GL 45	S* 60	PP	704348
GL 32	GL 25	ETFE	7043 75
GL 32	GL 28/ S 28	ETFE	7043 78
GL 32	GL 45	ETFE	7043 95
GL 45	GL 32	ETFE	7043 98
GL 45	GL 38	ETFE	7043 99
GL 45	S* 40	PTFE	7043 91
GL 32	NS 19/26	PP	704419
GL 32	NS 24/29	PP	704424
GL 32	NS 29/32	PP	704429

* buttress thread

Discharge tubes without recirculation valve

Discharge tube without recirculation valve for Dispensette® S

Packaging unit 1 pc.

Nominal volume ml	Model	Length (mm)	Order no.
1, 2, 5, 10	Fine tip	108	708002
5, 10	Standard	108	708005
25, 50, 100	Fine tip	135	708006
25, 50, 100	Standard	135	708008

Discharge tubes without recirculation valve for Dispensette® S Organic

Packaging unit 1 pc.

Nominal volume ml	Model	Length (mm)	Order no.
5, 10	Fine tip	108	708012
5, 10	Standard	108	708014
25, 50, 100	Fine tip	135	708016
25, 50, 100	Standard	135	708019

Discharge tubes with recirculation valve

Discharge tubes with recirculation valve for Dispensette® S

Packaging unit 1 pc.

Nominal volume ml	Model	Length (mm)	Order no.
1, 2, 5, 10	Fine tip	108	708102
5, 10	Standard	108	708104
25, 50, 100	Fine tip	135	708106
25, 50, 100	Standard	135	708109

Discharge tubes with recirculation valve for Dispensette® S Organic

Packaging unit 1 pc.

Nominal volume ml	Model	Length (mm)	Order no.
5, 10	Fine tip	108	708112
5, 10	Standard	108	708114
25, 50, 100	Fine tip	135	708116
25, 50, 100	Standard	135	708119

Flexible discharge tubing with recirculation valve for Dispensette® S and Dispensette® S Organic

PTFE, spiraled, approx. 800 mm long with safety handle.

Packaging unit 1 pc.

Not suitable for hydrofluoric acid (HF)

Nominal volume ml	Discharge tubing outer diameter	Discharge tubing inner diameter	Order no.
1, 2, 5, 10	3	2	708132
25, 50, 100	4.5	3	708134

Discharge valve Dispensette® S

PFA/Boro 3.3/ceramic/platinum-iridium.

No valve marking.

Packaging unit 1 pc.

for nominal volume ml	Order no.
1, 2*	6749
5, 10	6727
25, 50, 100	6728

* with valve marking '1 + 2'

Dispensette® S Organic discharge valve

PFA/Boro 3.3/ceramic/tantalum.

Valve marked with 'ORG'.

Packaging unit 1 pc.

for nominal volume ml	Order no.
5, 10	6729
25, 50, 100	6730

Dispensette® S and Dispensette® S Organic filling valve

PFA/ETFE/Boro 3.3/ceramic.


No valve marking.

Packaging unit 1 pc.

for nominal volume ml	Order no.
1, 2, 5, 10	6734

for nominal volume ml	Order no.
25, 50, 100	6735

Telescopic filling tubes for Dispensette® S and Dispensette® S Organic


FEP.

Individually adjustable lengths.

Packaging unit 1 pc.

for nominal volume ml	Outer diameter mm	Length (mm)	Order no.
1, 2, 5, 10	6	70–140	708210
		125–240	708212
		195–350	708214
		250–480	708216
25, 50, 100	7.6	170–330	708218
		250–480	708220

Venting stopper for microfilter with Luer cone

Description	Packaging unit	Order no.
Venting stopper for microfilter with Luer cone. PP. Venting stopper and PTFE-sealing ring.	1 pc.	704495

Recirculation tube

Description	Packaging unit	Order no.
Recirculation tube. FEP	1 pc.	6747

Adjustment/installation tool

Description	Packaging unit	Order no.
Adjustment/in-stallation tool	1 pc.	6748

Bottle stand


Description	Packaging unit	Order no.
Bottle stand. PP. Support rod, 325 mm, base plate 220 × 160 mm.	1 pc.	704275

Sealing ring for valve block

Description	Packaging unit	Order no.
Sealing ring for valve block. PTFE, for highly volatile media.	1 pc.	704486

Drying tube

Description	Packaging unit	Order no.
Drying tube, including PTFE sealing ring (without drying agent)	1 pc.	707930


Screw cap with fastener

Packaging unit 1 pc.

Description	Nominal volume ml	Order No.
PP, red, for Dispensette® S	1, 2, 5, 10	706018
	25, 50, 100	706019
PP, yellow, for Dispensette® S Organic	5, 10	706025
	25, 50, 100	706027
ETFE, for Dispensette® S and Dispensette® S Organic	1, 2, 5, 10	706029
PTFE, for Dispensette® S and Dispensette® S Organic	25, 50, 100	706031

Dispensette® drum extraction system

Packaging unit 1 pc.

Description	Order no.
Dispensette® drum extraction system	704261
Tripod connection for wall mount, stainless steel, Dispensette® extraction system	704268
Table/rack clamp for wall mount for Dispensette® extraction system	704272

15 Repairs

15.1 Sending for repair

NOTICE

Transporting of hazardous materials without a permit is a violation of federal law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Outside USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

16 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individ-

ual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The detailed testing instructions are available for download on www.brand.de or www.brandtech.com.

BRAND also offers the possibility to have your devices calibrated by our factory calibration service or by the BRAND DAkkS laboratory. Simply send us the device to be calibrated, accompanied by details about which type of calibration you would like. The device will be returned to you after a few days together with a test report (factory calibration) or a DAkkS calibration certificate. More information can be obtained from your dealer or directly from BRAND. The order document is available for download on www.brand.de (see Technical Documents).

17 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

18 Disposal

Before disposal, observe the relevant national disposal regulations, and ensure that the product is disposed of properly.

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