JULABO GmbH

Gerhard-Juchheim-Str. 1 Tel. +49 7823 51-0
77960 Seelbach Fax +49 7823 24 91
Germany info.de@julabo.com
www.julabo.com

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Kältebad-Grundgerät 1000FW

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Authorized representative in charge of administering technical documentation:

HR: Torsten Kauschke, im Haus / *on the manufacturer's premises as defined above*

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1000FW Basic Refrigerated Bath Units

Intended use

The basic refrigerated bath units can be combined with JULABO circulators from the CORIO CD, CORIO CP, DYNEO DD, and MAGIO MS product series. When combined with these circulators, they are designed to control the temperature of liquid media recommended by JULABO.

How to use the basic refrigerated bath unit is described in the circulator operating instructions.

⚠ WARNING

Danger from mains voltage!

Risk of injury from electric current.

- Compare the mains voltage and frequency with the details on the model plate.
- Only connect the unit to a fused mains connection via a FI circuit breaker ($I_a=30$ mA).
- The unit may only be connected to mains power outlets with a protective contact (PE)!
- The mains plug is a safe isolator from the power supply network and must be freely accessible at all times.
- Do not start the device if it has a damaged power cable.
- Check the mains cable regularly for signs of damage.
- We cannot accept any liability for incorrectly connecting to the power supply.



⚠ CAUTION

Cold or hot unit surfaces

Risk of frostbite, burns

- Device part surfaces can become very cold or hot
- Before draining the device, temper it to room temperature
- Wear personal protective equipment (gloves, safety shoes, goggles).

- 1 CAN bus cable for connecting to the circulator.
- 2 Recessed grip
- 3 Mains extension cable, cooling machine – circulator
- 4 CAN connector
- 5 Mains fuses
- 6 Mains connection for circulator power supply
- 7 Mains connection for cooling machine
- 8 Cooling water outlet, G3/4 AG
- 9 Cooling water inlet, G3/4 AG
- 10 Recessed grip
- 11 Removable ventilation grid (behind drain valve)

Included in delivery

- Mains cable and mains extension cable (N5/15)
- CAN bus connection cable

Technical data

Refrigeration basic units		1000FW					
Working temperature range	°C	-40 ... 200					
Cooling capacity (ethanol medium)	°C	20	0	-10	-20	-30	-40
	kW	1.00	0.98	0.75	0.53	0.27	0.13
Refrigerant		R449A					
Cooling water pressure inlet, max.	bar	6					
Cooling water differential pressure, min.	bar	0.3					
Cooling water consumption, min.	l/min	1.1					
Dimensions (WxDxH)	cm	42 x 49 x 70					
Useful bath opening (WxD)	cm	22 x 15					
Bath depth	cm	15					
Filling volume, from ... to	l	5.0 ... 7.5					
Weight	kg	48.5					
Protection class according to IEC 60529		IP21					
Ambient temperature range	°C	5 ... 40					
Mains connection	V/Hz	230 ±10% /50					
Power consumption (230 V)	A	Nom. 6 / Tot. 16					
CH model (at 230 V)	A	Nom. 6 / Tot. 10					
UK model (at 230 V)	A	Nom. 6 / Tot. 13					
Mains connection	V/Hz	208-230 ±10% /60					
Current consumption	A	Nom. 6 / Tot. 16					
Mains connection	V/Hz	115 ±10% / 60					
Current consumption	A	Nom. 9 / Tot. 16					
Mains connection	V/Hz	200 -5%; + 10% /50/60					
Current consumption	A	Nom. 11 / Tot. 15					

All data refers to rated voltage and rated frequency. Ambient temperature 20°C. Subject to technical changes.

Material of parts that come into contact with the medium

Part	Material
Bath tank and drain cock	1.4301/304H
Bath tank/bath cover gasket	FKM Viton®
O-ring, drain cock	FKM Viton®

Cooling water

Check the quality of the cooling water before connecting the device to a cooling water supply.

	NOTE
	<p>Corrosion and calcification caused by inadequate cooling water quality</p> <p>Inadequate cooling water quality can lead to calcification or corrosion of important system components and thus to failure of the device.</p> <ul style="list-style-type: none"> Observe the recommended cooling water quality Avoid hard water, water that contains iron or chlorine, distilled or de-ionized water, seawater or water that comes from cooling towers or rivers

Recommended cooling water quality:

Parameter	Permissible value
pH value	7.5 to 9.0
Sulfate [SO ₄ 2-]	< 100 ppm
Hydrogen carbonate [HCO ₃ ⁻] / sulfate [SO ₄ 2-]	> 1 ppm
Hardness [Ca ²⁺ , Mg ²⁺] / [HCO ₃ ⁻]	> 0.5 dH
Alkalinity	60 ppm < [HCO ₃ ⁻] < 300 ppm
Conductivity	< 500 µs / cm
Chloride (Cl ⁻)	< 50 ppm
Phosphate (PO ₄ ³⁻)	< 2 ppm
Ammoniac (NH ₃)	< 0.5 ppm
Free chlorine	< 0.5 ppm
3-valent iron ions (Fe ³⁺)	< 0.5 ppm
2-valent manganese ions (Mn ²⁺)	< 0.05 ppm
Carbon dioxide (CO ₂)	< 10 ppm
Hydrogen sulfide (H ₂ S)	< 50 ppm
Oxygen content	< 0.1 ppm
Algae growth	not permissible
Suspended solids	not permissible

- Recommended cooling water temperature range: 10°C ... 15°C
- Permissible cooling water temperature range: 5°C ... 30°C

If the recommended cooling water temperature range is exceeded, performance restrictions are to be expected. If the cooling water temperature is low, condensation may occur on the cooling circuit.

Connect cooling water to the device

The cooling water lines are mounted directly on the thread of the cooling water connections.

NOTE

Environmental damages due to leakage in the cooling water circuit

If there is a leak in the cooling water circuit, oil from the refrigerant cycle can get in to the cooling water and lead to environmental damages.

- Observe the legal regulations applicable at the place of operation
- If there is leakage in the cooling water circuit, shut down the device and repair the fault.

- ▶ The device is switched off.
 1. Connect the cooling water line to the cooling water inlet.
 2. Connect the cooling water return to the cooling water outlet.
 3. Route the line connected to the cooling-water outlet to the on-site drain or return flow.
- ✓ Cooling water is connected. After start-up, check all hose connections for leaks.

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