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Operation of manual crimping tools (Standard)

Crimpers from MACHEREY-NAGEL are supplied with a basic adjustment for the individual cap size and the septa thickness that is typically used for these caps. Crimp necks and septa thicknesses, however, can differ from manufacturer to manufacturer. Furthermore it needs especially to be considered with regard to 20 mm crimp caps that these are used in combination with vials having a flat DIN crimp neck as well as with those having a beveled top Headspace crimp neck (with a very low bearing surface for the septa).

Beveled Top Headspace Crimp Neck

Flat DIN Crimp Neck

In order to achieve an optimal crimp result, a further adjustment might be necessary. This may be conducted by the crimping pressure and/or the fixation of the crimping height.

An **optimal crimp result** is characterized by the following features:


1. The cap surface is flat and the cap sides fit tightly and firmly around the glass edges without any deformation of the aluminium surface.
2. The septa neither emerges out of the center hole of the cap nor is being sucked into the vial.
3. The cap cannot be turned by usage of single fingers without enormous expenditure of energy.

Often the so-called "twist test" of the crimp cap is being used, in order to verify the crimp result. This test is only conditionally meaningful. On one hand the slippery PTFE lamination of the septa results already by a relatively low expenditure of energy in a torque that enables the turn of the cap, even if the cap was crimped correctly. As the test often is done with the complete hand (instead of using single fingers) and with a high expenditure of energy, the false statement is even strengthened.

On the other hand overcrimped closures that cannot be turned, which have, however, either deformed cap sides or sucked in septa, may not be as tight as you would expect from such an "unturnable" crimp. This is due to the fact that septa material is sucked into the vial opening by the far too high crimping pressure and is then, however, missing at the cap sides for sealing. The same is also valid for septa that - due to a wrong adjustment of the crimping height and a too strong crimping pressure - emerge through the cap's center hole.

Construction of crimpers N 11, N 13 and N 20 (N 8 is only adjustable in crimping pressure)

1. Screw in the handle for the adjustment of the crimping pressure by the suspension point of the handles
2. Crimping head with a hexagon screw inside. With the attached hexagon key an adjustment of the crimping height may be carried out
3. Striker

Correction of inaccurate crimp features

- Crimp too loose (undercrimped) → Adjust crimping height (in the crimping head, see ②) and/or crimping pressure (by screw in the handle, see ①)
- Septa emerges through the center hole (overcrimped) → Adjust crimping height (in the crimping head, see ②) and/or crimping pressure (by screw in the handle, see ①)
- Deformation of the cap sides (overcrimped) → Adjust crimping height (in the crimping head, see ②)
- Septa is sucked into the vial (overcrimped) → Adjust crimping pressure (by screw in the handle, see ①)

Adjustment of the crimping pressure

1. Further unscrew the screw ① in the handle → handles cannot be pressed so close together → lower crimping pressure
2. Further screw the screw ① in the handle → handles can be pressed closer together → higher crimping pressure

Adjustment of the crimping height

1. Fixate the handles of the crimper with the striker ③.
2. Hold crimping head ② firmly in one hand and insert the supplied hexagon key in the hexagon screw inside the crimping head.
3. a) Make a half up to a complete turn **to the left (counter clockwise rotation), in case the crimp is too loose.**
 b) Turn the hexagon key **to the right (clockwise rotation), in case the crimp is too tight.**

Ordering information	REF
Manual crimpers, Standard (adjustable in crimping height and crimping pressure)	
Crimper for 8 mm aluminium crimp caps (only pressure adjustable)	735126
Crimper, height adjustable, for 11 mm crimp caps	735111
Crimper, height adjustable, for 13 mm crimp caps	735113
Crimper, height adjustable, for 13 mm Flip Top / Flip Off caps	735133
Crimper, height adjustable, for 20 mm crimp caps	735120
Crimper, height adjustable, for 20 mm Flip Top / Flip Off caps	735132
Manual decappers, Standard (do not require any adjustment)	
Decapper for 8 mm crimp caps	735408
Decapper for 11 mm crimp caps	735911
Decapper for 13 mm crimp caps	735913
Decapper for 20 mm crimp caps	735920

For further manual/ electronic crimping tools please see our website under www.mn-net.com.