

# NucleoSpin® Plasmid (NoLid)

## 1 Kit contents

NucleoSpin® Plasmid (NoLid)		
REF	12 preps 740499.12C	240 preps 740499.240C
Resuspension Buffer A1	15 mL (add RNase A before first use)	100 mL (add RNase A before first use)
Lysis Buffer A2	15 mL	100 mL
Neutralization Buffer A3	20 mL	200 mL
Wash Buffer AW	30 mL	2 x 100 mL
Wash Buffer A4 (Concentrate)	6 mL (add 24 mL ethanol before first use)	100 mL (add 400 mL ethanol before first use)
Elution Buffer AE	30 mL	125 mL
RNase A (lyophilized)	6 mg	40 mg
NucleoSpin® Plasmid (NoLid) Columns (white rings)	12	4 x 60
User manual	1	1

## 2 How to use the kit

Please see the protocol information how to use the kit (see next pages). For further questions and more detailed information, please contact MACHEREY-NAGEL at [tech-bio@mn-net.com](mailto:tech-bio@mn-net.com) for protocol information how to use the kit on specific robotic instruments.

For storage conditions, product use restrictions, and safety information, please see the general NucleoSpin® Plasmid (NoLid) user manual.

### NucleoSpin® kits on QIAcube®

MN is not recommending to use this kit on specific robots. The use of NucleoSpin® kits on the QIAcube® is solely at your own discretion. MACHEREY-NAGEL is not responsible for loss of warranty claims or other consequences.

### 3 General information

Application:	Plasmid DNA
Kit:	NucleoSpin® Plasmid (NoLid) (REF 740499.240C) instead of: QIAprep® Spin Miniprep Kit
Sample material:	Up to 5 mL LB cultures (for TB or 2x YT medium use up to 2.5 mL culture volume)
Protocol name:	Plasmid DNA purification using the QIAprep® Spin Miniprep Kit (with PB wash) DNA_QIAprepMiniprep_UpTo5mLLBCultures_Standard_V3
Editable parameters:	Elution volume: 30–100 µL; default 50 µL

### 4 Using the kit


1. Fill the designated buffer bottles with the buffers according to the buffer table below (make sure that RNase A was added to Buffer A1 in advance).
2. Harvest bacterial cells of up to 5 mL LB cultures by centrifugation at 6800 x g for 3 minutes.
3. Place Sample Tubes (2 mL Safe-Lock microcentrifuge tubes) containing the pelleted bacterial cells into sample rack (shaker).
4. Insert disposable Filter Tips 1000 µL and 1000 µL wide-bore.
5. General equipment setup is shown below.

### 5 Additional materials

Refer to the QIAcube® protocol sheet for required consumables (e.g., sample tubes, collection tubes, instrument accessories, disposable tips, etc.) and software requirements.

## 6 Rotor adapter

Position	Labware	Lid position
1	NucleoSpin® Plasmid (NoLid) Column	L1
2	–	–
3	1.5 mL collection tube*	L3



\* Sarstedt, Micro tube 1.5 mL Safety Cap

## 7 Buffers (Reagent Bottle Rack)

Position	MN Reagent	Replaced QIAGEN® Reagent
1	Buffer A1	Buffer P1
2	Buffer A2	Buffer P2
3	Buffer A3	Buffer N3
4	Buffer AW	Buffer PB
5	Buffer A4	Buffer PE
6	Buffer AE	Buffer EB

### NucleoSpin® kits on QIAcube®

MN is not recommending to use this kit on specific robots. The use of NucleoSpin® kits on the QIAcube® is solely at your own discretion. MACHEREY-NAGEL is not responsible for loss of warranty claims or other consequences.



## 8 Disclaimer

MACHEREY-NAGEL GmbH & Co. KG makes every effort to include accurate and up-to-date information within this publication; however, it is possible that omissions or errors might have occurred. MACHEREY-NAGEL GmbH & Co. KG cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided in this publication. Changes in this publication can be made at any time without notice. For technical details and detailed procedures of the specifications provided in this document please contact your MACHEREY-NAGEL representative. This publication may contain reference to applications and products which are not available in all markets. Please check with your local sales representative.

All mentioned trademarks are protected by law. All used names and denotations can be brands, trademarks, or registered labels of their respective owner – also if they are not special denotation. To mention products and brands is only a kind of information (i.e., it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment). Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or support by MACHEREY-NAGEL GmbH & Co. KG. Any views or opinions expressed herein by the authors' do not necessarily state or reflect those of MACHEREY-NAGEL, and shall not be used for advertising or product endorsement purposes.

Trademarks: QIAGEN®, QIAcube®, and QIAamp® are registered trademarks of the QIAGEN Group;

NucleoSpin is a registered trademark of MACHEREY-NAGEL GmbH & Co. KG.

