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
## Operation of manual, ergonomic crimpers

Congratulations to the purchase of a manual crimping tool of the latest generation that especially has been optimized with regard to ergonomics and handling. In comparison to the standard crimping tools these crimpers have the following advantages:

- Regulation of the crimp set by the adjustment knob on top of the crimping head. The knob is easily accessible and visible.
- The conducted adjustment doesn't need any further fixation (e.g. by a fixing nut)
- The ergonomic crimping tools are much lighter than the all-steel crimpers, however, with the same level of robustness, as the body is constructed of tough, light glass filled plastic and with integrated handle reinforcements.
- The curved, ergonomically shaped handles are more comfortable in the hand in comparison to the block metal handles of the standard crimping tools.
- Bottom-pull handle makes it easier to hold the head steady during crimping in comparison to the standard crimpers that require pushing the handles together for crimping.
- Due to the long, rectangular arranged crimping head you have a better vertical clearance over the vial.
- The stroke of the crimper can be seen by watching the adjustment knob that touches the crimper body once the crimp setting has been reached and the crimping process stops (you don't have to keep squeezing after the stop position has been reached).
- A label on the front of the crimper (not only on the packaging) allows identification of brand and product.

Crimpers from MACHEREY-NAGEL are supplied with a basic adjustment for the individual cap size and the septa thickness that is typically used for these caps. Crimp necks and septa thicknesses, however, can differ from manufacturer to manufacturer. Furthermore it needs especially to be considered with regard to 20 mm crimp caps that these are used in combination with vials having a flat DIN crimp neck as well as with those having a beveled top Headspace crimp neck (with a very low bearing surface for the septa).

In order to achieve an optimal crimp result, a further adjustment might be necessary. With the ergonomic crimpers this is carried out by the adjustment knob on top of the crimping head.

The adjustment knob, which sets the stroke of the tool and the tightness of the crimp, is conveniently positioned at the head of the crimper. The easily-viewed knob clearly displays + and - symbols along with directional arrows, in order to simplify adjusting the level of crimp desired. During the crimping process, in which the bottom handle is pulled up, the adjustment knob moves towards the crimper body. Once the knob touches the head, the crimp setting has been reached and the crimping process stops.

An **optimal crimp result** is characterized by the following features:

1. The cap surface is flat and the cap sides fit tightly and firmly around the glass edges without any deformation of the aluminium surface.
2. The septa neither emerges out of the center hole of the cap nor is being sucked into the vial.
3. The cap cannot be turned by usage of single fingers without enormous expenditure of energy.

Often the so-called „twist test“ of the crimp cap is being used, in order to verify the crimp result. This test is only of limited relevance. On one hand the slippery PTFE lamination of the septa results already by a relatively low expenditure of energy in a torque that enables the turn of the cap, even if the cap was crimped correctly. As the test often is done with the complete hand (instead of using single fingers) and with a high expenditure of energy, the false statement is even strengthened.

On the other hand overcrimped closures that cannot be turned, which have, however, either deformed cap sides or sucked in septa, may not be as tight as you would expect from an optimal crimp. This is due to the fact that septa material is sucked into the vial opening by the far too high crimping pressure and is then, however, missing at the cap sides for sealing. The same is also valid for septa that – due to a wrong adjustment of the crimping height and a too strong crimping pressure – emerge through the cap's center hole.

**Please consider that manual decappers don't need any adjustment. The adjustment knob on top of the decapper has no function at all.**

Ordering information	REF
<b>Manual crimpers</b> (ergonomic, crimping pressure adjustable by knob on the crimping head)	
Manual, ergonomic crimper for 8 mm crimp caps	735208
Manual, ergonomic crimper for 11 mm crimp caps	735211
Manual, ergonomic crimper for 20 mm crimp caps	735220
<b>Manual decappers</b> (ergonomic, don't require any adjustment despite the existing adjustment knob)	
Manual, ergonomic decapper for 11 mm crimp caps	735311
Manual, ergonomic decapper for 20 mm crimp caps	735320
<b>Accessories for ergonomic crimping tools</b>	
Rack for two ergonomic crimping tools, blue, 240 x 95 x 65 mm	735509

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**MACHEREY-NAGEL GmbH & Co. KG**
**MACHEREY-NAGEL AG**