

Now sold under the
Thermo Scientific brand

Thermo
SCIENTIFIC
Part of Thermo Fisher Scientific

Acclaim Explosives E1 and E2 Columns Quick Start

Acclaim Explosives E1, 5 µm, Analytical Column, 4.6 x 250 mm (P/N 064305)
 Acclaim Explosives E2, 5 µm, Analytical Column, 4.6 x 250 mm (P/N 064309)
 Acclaim Explosives E2, 3µm, Analytical Column, 3.0 x 250 mm (P/N 070081)
 Acclaim Explosives E2, 3µm, Analytical Column, 3.0 x 150 mm (P/N 070082)
 Acclaim Explosives E2, 3µm, Analytical Column, 2.1 x 150 mm (P/N 070083)
 Acclaim Explosives E2 RSLC, 2.2µm, Analytical Column, 2.1 x 100mm (P/N 076225)
 Acclaim Explosives E2 RSLC, 2.2µm, Analytical Column, 2.1 x 150 mm (P/N 076226)
 Acclaim Explosives E2 RSLC, 2.2µm, Analytical Column, 3.0 x 100 mm (P/N 076227)
 Acclaim Explosives Kit (E1 and E2 5 µm, 4.6 x 250 mm Analytical Columns and Guards) (P/N 064312)

1. Overview

The Acclaim® Explosives E1 and E2 columns are silica-based reversed phase columns designed specifically for separating the 14 explosives compounds listed in EPA SW-846 Method 8330. Both columns exhibit superior resolution and robustness compared with other commercially available columns, with selectivity that is complementary to each other.

2. Main features

- Baseline resolution of all 14 compounds targeted by EPA Method 8330
- Mutually-complementary selectivity
- Simple, isocratic elution conditions
- Rugged columns with good lot-to-lot reproducibility
- Unique selectivities for separating other nitroaromatic molecules

3. Physical data

Column Chemistry: Proprietary
 Substrate: High-purity, porous, spherical silica
 Particle size: 2.2, 3, and 5 µm
 Surface area: 300 m²/g
 Pore size: 120 Å

4. Specifications and Recommended Operational Parameters

Shipping solution: 70% acetonitrile in water
 Long term storage solution: 100% acetonitrile or methanol
 Aqueous compatibility: 100% aqueous compatibility

Column Type	Particle Size	Column Dimensions	Maximum Pressure (bar)	Maximum Flow Rate (mL/min)	Typical Flow Rate (mL/min)	pH Range	Typical Temperature (°C)	Maximum Temperature (°C)
Acclaim E1	5µm	4.6 x 250mm	400	2.0	0.8 – 1.5	3.5 – 7.0	25 – 35	50
Acclaim E2	5µm	4.6 x 250mm	400	2.0	0.8 – 1.5	2.5 – 7.5	25 – 35	50
Acclaim E2	3µm	3.0 x 250mm	800	1.0	0.4 – 0.6	2.5 – 7.5	25 – 35	50
		3.0 x 150mm	600	1.3	0.4 – 0.8	2.5 – 7.5	25 – 35	50
		2.1 x 150mm	600	0.8	0.2 – 0.5	2.5 – 7.5	25 – 35	50
Acclaim E2	2.2µm	3.0 x 100mm	600	1.6	0.4 – 1.4	2.5 – 7.5	25 – 35	50
		2.1 x 150mm	800	1.0	0.25 – 0.7	2.5 – 7.5	25 – 35	50
		2.1 x 100mm	700	1.0	0.25 – 0.7	2.5 – 7.5	25 – 35	50

5. Operational Guidelines

- The direction of flow is marked on the column. While it is not harmful to the column, reverse flow should be avoided except to attempt removal of inlet blockage (see “Column Care” in the Manual).
- A new column is shipped in a solution of acetonitrile and water in 70:30 ratio. Initially, care should be taken to avoid any mobile phase that might cause a precipitate.
- The use of a guard column and/or a pre-column filter is highly recommended whenever possible to protect the analytical column and extend its useful lifetime.
- Recommended operating pressures are column format dependent (see Table above for details).
- Columns must not be exposed to a condition outside the recommended pH range or temperature range; permanent changes in selectivity may occur in such cases.

Suggested conditions for EPA Method 8330. Refer to the manual for further detail and advice on optimization.

Column Type	Particle Size	Column Dimensions	Mobile phase	Temperature	Flow Rate
Acclaim E1	5 µm	4.6 x 250 mm	43% methanol	31 °C	1.0 mL/min
Acclaim E2	5 µm	4.6 x 250 mm	48% methanol	30 °C	1.0 mL/min
Acclaim E2	3 µm	3.0 x 250 mm	49% methanol	27 °C	0.54 mL/min
		3.0 x 150 mm	48% methanol	30 °C	0.54 mL/min
		2.1 x 150 mm	48% methanol	30 °C	0.26 mL/min
Acclaim E2	2.2 µm	3.0 x 100 mm	48% methanol	31 °C	0.70 mL/min
		2.1 x 150 mm	48% methanol	31 °C	0.34 mL/min
		2.1 x 100 mm	48% methanol	31 °C	0.34 mL/min

Other applications for these columns may be found in the Manual or in the Acclaim Library.