Optimisation of solid phase extraction for the analysis of benzodiazepines from plasma

¹Y. T. Hui, ²T. Liddicoat, ²L. Pereira, ²T. Edge

¹Department of Forensic Science and Drug Monitoring, Kings Collegue London, UK,

²Thermo Fisher Scientific, Runcorn, Cheshire, UK

Introduction

Benzodiazepines possess sedative, hypnotic, anxiolytic, anticonvulsant, muscle relaxant, and amnesic actions, which are useful in a variety of indications such as alcohol dependence, seizures, anxiety, panic, agitation and insomnia. This class of compounds has also been associated with drug abuse.

The analysis of benzodiazepines within biological systems allows effective treatment of a range of anxiety based conditions. This requires a combination of mass spectrometry, chromatography and also sample preparation to allow the quantitation of these compounds within the biological fluid, typically plasma. The sample preparation has a significant effect on the quality of the data being produced. The subsequent sections detail the investigation into various sample preparation techniques for the removal of endogenous matrix components from the original plasma sample.

Methods

Annroaches

Two sample preparation techniques were investigated with the aim of determining the matrix components left after the extraction process:

1. Protein Precipitation

100 μ l of plasma was added to 300 μ l of ACN. Supernatant blown down and reconstituted in 500 μ L of water-acetonitrile-formic acid (95:5:0.1 ν / ν)

2. Solid Phase extraction

Details of the optimisation procedure are given below using either the HyperSEP Retain PEP, or the HyperSEP Retain CX solid phase extraction sorbents. Optimisation of the washes is given in the subsequent section

Chromatographic and Detection Methods

Samples were prepared using the two different sample preparation techniques and then analysed using LC-MS. Instrumentation: Thermo Scientific Accela UHPLC with Thermo Scientific MSQ mass spectrometer.

Column – Hypersil GOLD 1.9µm, 50 x 2.1 mm

Mobile Phase A - 0.1% Formic Acid (aq); Mobile Phase B - 0.1% Formic Acid in ACN.

Gradient: 95% [A] and 5% [B] for 0.5 min, linear change to 5% [A] and 95% [B] in 3 min, then change to 95% [A] and 5% [B] and hold for 1 min.


Flow rate: 0.6mL/min; Injection volume: 10µL


Detection: positive ESI; for elution profile, 12 benzodiazepines SIM for each $[M+H]^{*}$; for blank plasma extracts, full scan from m/z 100 - 1000

Elution Profiles

This offers a quick optimisation of wash and elution conditions and recovery.


Compounds are loaded onto a single SPE cartridge and washed with increasing elutropic solvent strength (Table 1). Flow injection analysis of the eluents will determine the optimal extraction conditions, Figure 1.





Results

Figure 2 shows the experimental data for the optimisation of the benzodiazepines extraction on the HyperSEP Retain CX and HyperSEP Retain PEP materials. From this graph the optimum wash and elution conditions for all of the compounds can be readily determined.





Matrix components after protein precipitation

Analysis of the matrix was achieved using TIC chromatograms. Since the chromatography was not optimised for the matrix components, several aqueous blank injections cycles were also performed subsequent to the plasma extract injection to elute highly retained components. This is important information as it is generally assumed that the column has not altered from one injection to the next.

assumed that the column has not altered from one injection to the next.


It can be seen from Figure 3 that there are substantial amounts of matrix components still remaining in the sample after protein precipitation, in particular there are some very defined regions, of early eluting polar compounds, and then more hydrophobic compounds which elute at later retention times.





Optimising the wash conditions

The use of generic methods in SPE with a 5 or 10% organic wash results in a substantial amount of matrix components still remaining in the injected samples. Increasing the wash concentration to 60% methanol removed the majority of the polar compounds but did not remove the hydrophobic components, which required up to 80 or 90% methanol washing before the extract became clean.



Conclusions

The importance of the sample preparation technique in removing endogenous matrix components has been demonstrated. It has also been shown that with typical methods used in routine bioanalysis it takes several injection cycles to remove the endogenous matrix components from the system. This suggest that the chromatographic system changes on the first few injections before stabilising.

A quick approach to the optimisation of solid phase extraction method has been presented with two proposed optimised methods for the determination of a series of twelve benzodiazepines.

For additional information, please visit our Chromatography Resource Centre which can be found at:

©2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries .

