

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 02/19/2020

Reviewed on 02/19/2020

**1 Identification**

- **Product identifier**
- **Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**
- **Part number:** 5190-9767
- **Application of the substance / the mixture** Reagents and Standards for Analytical Chemical Laboratory Use
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Agilent Technologies, Inc. Tel: 800-227-9770  
5301 Stevens Creek Blvd  
Santa Clara, CA 95051,  
USA
- **Information department:** e-mail: [pdl-msds\\_author@agilent.com](mailto:pdl-msds_author@agilent.com)
- **Emergency telephone number:** CHEMTREC®: 1-800-424-9300

**2 Hazard(s) identification**

- **Classification of the substance or mixture**



GHS05 Corrosion

Met. Corr.1 H290 May be corrosive to metals.



GHS07

Skin Irrit. 2 H315 Causes skin irritation.

Eye Irrit. 2A H319 Causes serious eye irritation.

- **Label elements**

- **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

- **Hazard pictograms**



GHS05

- **Signal word** Warning

- **Hazard statements**

H290 May be corrosive to metals.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

- **Precautionary statements**

P280 Wear protective gloves / eye protection / face protection.

P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P321 Specific treatment (see on this label).

P332+P313 If skin irritation occurs: Get medical advice/attention.

(Contd. on page 2)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 1)

P337+P313 *If eye irritation persists: Get medical advice/attention.*  
 P406 *Store in corrosive resistant container with a resistant inner liner.*

- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**

HEALTH	2	Health = 2
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

### 3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:**  
Aqueous solution.  
Also contains substances at levels not considered to be hazardous.

- **Dangerous components:**

CAS: 7697-37-2	Nitric acid	<2%
RTECS: QU5775000	☠ Ox. Liq. 2, H272; ☠ Met. Corr.1, H290; Skin Corr. 1A, H314	

- **Additional information:**

The concentration of the acid stated in this SDS is calculated as an absolute mass concentration (%w/v). This is less than the acid concentration stated on the product label and COA, which reflects a percent value of the commercially available concentrated aqueous form of the acid.

### 4 First-aid measures

- **Description of first aid measures**
- **General information:** Immediately remove any clothing soiled by the product.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:**  
Immediately wash with water and soap and rinse thoroughly.  
If skin irritation continues, consult a doctor.
- **After eye contact:**  
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** Rinse mouth. Do not induce vomiting.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.

(Contd. on page 3)

## Safety Data Sheet

*acc. to OSHA HCS*

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 2)

- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.
- **Special hazards arising from the substance or mixture**  
Formation of toxic gases is possible during heating or in case of fire.
- **Advice for firefighters**
- **Protective equipment:** Wear self-contained respiratory protective device.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**  
Wear protective equipment. Keep unprotected persons away.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**  
Use neutralizing agent.  
Dispose contaminated material as waste according to item 13.  
Absorb liquid components with liquid-binding material.  
**DO NOT USE SAWDUST.**
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

- **PAC-1:**

CAS: 7697-37-2	Nitric acid	0.16 ppm
CAS: 7440-02-0	Nickel	4.5 mg/m <sup>3</sup>
CAS: 7440-22-4	Silver, powder	0.3 mg/m <sup>3</sup>
CAS: 7782-49-2	Selenium	0.6 mg/m <sup>3</sup>
CAS: 1314-62-1	vanadium pentoxide	0.64 mg/m <sup>3</sup>
CAS: 7440-48-4	Cobalt	0.18 mg/m <sup>3</sup>
CAS: 7440-28-0	thallium	0.06 mg/m <sup>3</sup>

- **PAC-2:**

CAS: 7697-37-2	Nitric acid	24 ppm
CAS: 7440-02-0	Nickel	50 mg/m <sup>3</sup>
CAS: 7440-22-4	Silver, powder	170 mg/m <sup>3</sup>
CAS: 7782-49-2	Selenium	6.6 mg/m <sup>3</sup>
CAS: 1314-62-1	vanadium pentoxide	7 mg/m <sup>3</sup>
CAS: 7440-48-4	Cobalt	2 mg/m <sup>3</sup>
CAS: 7440-28-0	thallium	3.3 mg/m <sup>3</sup>

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 02/19/2020

Reviewed on 02/19/2020

Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO3 [100ml bottle]

(Contd. of page 3)

Table with 3 columns: CAS number, Name, and Concentration. Rows include Nitric acid (92 ppm), Nickel (99 mg/m³), Silver, powder (990 mg/m³), Selenium (40 mg/m³), vanadium pentoxide (70 mg/m³), Cobalt (20 mg/m³), and thallium (20 mg/m³).

7 Handling and storage

- Handling:
Precautions for safe handling Store in cool, dry place in tightly closed receptacles.
Information about protection against explosions and fires: No special measures required.
Conditions for safe storage, including any incompatibilities
Storage:
Requirements to be met by storerooms and receptacles: Please refer to the manufacturers certificate for specific storage and transport temperature conditions.
Information about storage in one common storage facility: Store away from foodstuffs.
Further information about storage conditions: Keep receptacle tightly sealed.
Specific end use(s) No further relevant information available.

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters

Components with limit values that require monitoring at the workplace:

CAS: 7697-37-2 Nitric acid

Table with 2 columns: Exposure Limit and Value. Rows include PEL (Long-term value: 5 mg/m³, 2 ppm), REL (Short-term value: 10 mg/m³, 4 ppm; Long-term value: 5 mg/m³, 2 ppm), and TLV (Short-term value: 10 mg/m³, 4 ppm; Long-term value: 5.2 mg/m³, 2 ppm).

- Additional information: The lists that were valid during the creation were used as basis.
Exposure controls
Personal protective equipment:
General protective and hygienic measures: Keep away from foodstuffs, beverages and feed. Immediately remove all soiled and contaminated clothing. Wash hands before breaks and at the end of work. Avoid contact with the eyes and skin.

(Contd. on page 5)

## Safety Data Sheet

*acc. to OSHA HCS*

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 4)

· **Breathing equipment:**

Not required.

Use suitable respiratory protective device in case of insufficient ventilation.

· **Protection of hands:**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

The protective gloves to be used must comply with the specifications of EC Directive 89/686/EEC and the related standard EN374



Protective gloves

· **Material of gloves**

PVC gloves

Neoprene gloves

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:**



Tightly sealed goggles

### 9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

<b>Form:</b>	Liquid
<b>Color:</b>	Colorless
<b>Odor:</b>	Odorless
<b>Odor threshold:</b>	Not determined.

· **pH-value:** < 2

· **Change in condition**

<b>Melting point/Melting range:</b>	Not determined.
<b>Boiling point/Boiling range:</b>	Not determined. 100 °C (212 °F)

· **Flash point:** Not applicable.

· **Flammability (solid, gaseous):** Not determined.

· **Ignition temperature:** Not determined

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Product is not selfigniting.

(Contd. on page 6)

## Safety Data Sheet

*acc. to OSHA HCS*

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 5)

<b>· Danger of explosion:</b>	<i>Not determined.</i>
<b>· Explosion limits:</b>	
<b>Lower:</b>	<i>Not determined.</i>
<b>Upper:</b>	<i>Not determined.</i>
<b>· Vapor pressure at 20 °C (68 °F):</b>	<i>23 hPa (17.3 mm Hg)</i>
<b>· Density:</b>	<i>Not determined.</i>
<b>· Relative density</b>	<i>Not determined.</i>
<b>· Vapor density</b>	<i>Not determined.</i>
<b>· Evaporation rate</b>	<i>Not determined.</i>
<b>· Solubility in / Miscibility with Water:</b>	<i>Fully miscible.</i>
<b>· Partition coefficient (n-octanol/water):</b>	<i>Not determined.</i>
<b>· Viscosity:</b>	
<b>Dynamic:</b>	<i>Not determined.</i>
<b>Kinematic:</b>	<i>Not determined.</i>
<b>· Solvent content:</b>	
<b>VOC content:</b>	<i>0.00 %</i>
<b>· Other information</b>	<i>No further relevant information available.</i>

### 10 Stability and reactivity

- **Reactivity**  
*Stable under normal conditions.*  
*No further relevant information available.*
- **Chemical stability** *Stable under normal conditions.*
- **Thermal decomposition / conditions to be avoided:**  
*Formation of toxic gases is possible during heating or in case of fire.*
- **Possibility of hazardous reactions** *No dangerous reactions known.*
- **Conditions to avoid** *Heat.*
- **Incompatible materials:**  
*Strong oxidizing agents.*  
*Metals.*
- **Hazardous decomposition products:** *Formation of toxic gases is possible during heating or in case of fire.*

### 11 Toxicological information

- **Information on toxicological effects**
  - **Acute toxicity:**
- |   |   |
|---|---|
| <b>· LD/LC50 values that are relevant for classification:</b> |   |
| <b>CAS: 7440-22-4 Silver, powder</b>                          |   |
| <i>Oral</i>   | <b>LD50</b> <i>&gt;10,000 mg/kg (mouse)</i> |
- **Primary irritant effect:**
  - **on the skin:** *Irritant to skin and mucous membranes.*

(Contd. on page 7)

## Safety Data Sheet

*acc. to OSHA HCS*

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 6)

- **on the eye:** Irritating effect.

- **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:

Irritant

- **Carcinogenic categories**

- **IARC (International Agency for Research on Cancer)**

CAS: 7440-02-0	Nickel	2B
CAS: 7782-49-2	Selenium	3
CAS: 1314-62-1	vanadium pentoxide	2B
CAS: 7440-48-4	Cobalt	2B

- **NTP (National Toxicology Program)**

CAS: 7440-02-0	Nickel	R
CAS: 7440-48-4	Cobalt	R

- **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**

- **Aquatic toxicity:**

**CAS: 7440-22-4 Silver, powder**

LC50/24	0.015 mg/L (crustacean)
EC50/72h	0.00198 mg/l (Algae)
LC50/96 h	0.00807 mg/l (fish)

- **Persistence and degradability** No further relevant information available.

- **Behavior in environmental systems:**

- **Bioaccumulative potential** No further relevant information available.

- **Mobility in soil** No further relevant information available.

- **Additional ecological information:**

- **General notes:**

Not hazardous for water.

Must not reach bodies of water or drainage ditch undiluted or unneutralized.

- **Results of PBT and vPvB assessment**

- **PBT:** Not applicable.

- **vPvB:** Not applicable.

- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**

- **Recommendation:**

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

(Contd. on page 8)

## Safety Data Sheet

*acc. to OSHA HCS*

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 7)

- **Uncleaned packagings:**
- **Recommendation:** *Dispose in accordance with national regulations.*
- **Recommended cleansing agent:** *Water, if necessary with cleansing agents.*

### 14 Transport information

<ul style="list-style-type: none"> <li>· <b>UN-Number</b></li> <li>· <b>DOT, ADR, IMDG, IATA</b></li> <li>· <b>DOT</b></li> <li>· <b>ADR</b></li> <li>· <b>IMDG, IATA</b></li> </ul>	<p>UN3264</p> <p><i>Corrosive liquid, acidic, inorganic, n.o.s. (Nitric acid)</i></p> <p>3264 CORROSIVE LIQUID, ACIDIC, INORGANIC, N.O.S. (NITRIC ACID)</p> <p>CORROSIVE LIQUID, ACIDIC, INORGANIC, N.O.S. (NITRIC ACID)</p>
--	--

- **Transport hazard class(es)**

- **DOT**



<ul style="list-style-type: none"> <li>· <b>Class</b></li> <li>· <b>Label</b></li> </ul>	<p>8 Corrosive substances</p> <p>8</p>
--	--

- **ADR, IMDG, IATA**



<ul style="list-style-type: none"> <li>· <b>Class</b></li> <li>· <b>Label</b></li> </ul>	<p>8 Corrosive substances</p> <p>8</p>
--	--

- **Packing group**

<ul style="list-style-type: none"> <li>· <b>DOT, ADR, IMDG, IATA</b></li> </ul>	<p>III</p>
---	------------

<ul style="list-style-type: none"> <li>· <b>Environmental hazards:</b></li> </ul>	<p><i>Not applicable.</i></p>
---	-------------------------------

<ul style="list-style-type: none"> <li>· <b>Special precautions for user</b></li> </ul>	<p><i>Warning: Corrosive substances</i></p>
---	---

<ul style="list-style-type: none"> <li>· <b>Danger code (Kemler):</b></li> </ul>	<p>80</p>
--	-----------

<ul style="list-style-type: none"> <li>· <b>EMS Number:</b></li> </ul>	<p>F-A,S-B</p>
--	----------------

<ul style="list-style-type: none"> <li>· <b>Segregation groups</b></li> </ul>	<p>Acids</p>
---	--------------

<ul style="list-style-type: none"> <li>· <b>Stowage Category</b></li> </ul>	<p>A</p>
---	----------

<ul style="list-style-type: none"> <li>· <b>Stowage Code</b></li> </ul>	<p>SW2 <i>Clear of living quarters.</i></p>
---	---

<ul style="list-style-type: none"> <li>· <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b></li> </ul>	<p><i>Not applicable.</i></p>
--	-------------------------------

(Contd. on page 9)



## Safety Data Sheet

*acc. to OSHA HCS*

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 8)

**· Transport/Additional information:**

**· ADR**

**· Excepted quantities (EQ)**

Code: E1

Maximum net quantity per inner packaging: 30 ml

Maximum net quantity per outer packaging: 1000 ml

**· UN "Model Regulation":**

UN 3264 CORROSIVE LIQUID, ACIDIC, INORGANIC, N.O.S. (NITRIC ACID), 8, III

### 15 Regulatory information

**· Safety, health and environmental regulations/legislation specific for the substance or mixture**

**· Sara**

**· Section 355 (extremely hazardous substances):**

CAS: 7697-37-2	Nitric acid
CAS: 1314-62-1	vanadium pentoxide

**· Section 313 (Specific toxic chemical listings):**

CAS: 7697-37-2	Nitric acid
CAS: 7440-02-0	Nickel
CAS: 7440-22-4	Silver, powder
CAS: 7782-49-2	Selenium
CAS: 1314-62-1	vanadium pentoxide
CAS: 7440-48-4	Cobalt
CAS: 7440-28-0	thallium

**· TSCA (Toxic Substances Control Act):**

All components have the value ACTIVE.

**· Hazardous Air Pollutants**

CAS: 7440-48-4	Cobalt
----------------	--------

**· Proposition 65**

**· Chemicals known to cause cancer:**

CAS: 7440-02-0	Nickel
CAS: 1314-62-1	vanadium pentoxide
CAS: 7440-48-4	Cobalt

**· Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

**· Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

**· Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

(Contd. on page 10)

## Safety Data Sheet

*acc. to OSHA HCS*

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 9)

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

CAS: 7440-22-4	Silver, powder	D
CAS: 7782-49-2	Selenium	D

· **TLV (Threshold Limit Value established by ACGIH)**

CAS: 7440-02-0	Nickel	A5
CAS: 1314-62-1	vanadium pentoxide	A3
CAS: 7440-48-4	Cobalt	A3

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

CAS: 7440-02-0	Nickel	
----------------	--------	--

· **Hazard pictograms**



GHS05

· **Signal word** Warning

· **Hazard statements**

H290 May be corrosive to metals.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

· **Precautionary statements**

P280 Wear protective gloves / eye protection / face protection.

P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P321 Specific treatment (see on this label).

P332+P313 If skin irritation occurs: Get medical advice/attention.

P337+P313 If eye irritation persists: Get medical advice/attention.

P406 Store in corrosive resistant container with a resistant inner liner.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### 16 Other information

The information contained in this document is based on Agilent's state of knowledge at the time of preparation. No warranty as to its accurateness, completeness or suitability for a particular purpose is expressed or implied.

· **Date of preparation / last revision** 02/19/2020 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

(Contd. on page 11)

**Safety Data Sheet**  
*acc. to OSHA HCS*

Printing date 02/19/2020

Reviewed on 02/19/2020

**Product name: ICH/USP Oral Target Elements Standard B: Ni @ 200; Ag, Se @ 150; V @ 100; Co @ 50; Tl @ 8 µg/mL in 2% HNO<sub>3</sub> [100ml bottle]**

(Contd. of page 10)

*HMIS: Hazardous Materials Identification System (USA)*

*VOC: Volatile Organic Compounds (USA, EU)*

*LC50: Lethal concentration, 50 percent*

*LD50: Lethal dose, 50 percent*

*PBT: Persistent, Bioaccumulative and Toxic*

*vPvB: very Persistent and very Bioaccumulative*

*NIOSH: National Institute for Occupational Safety*

*OSHA: Occupational Safety & Health*

*TLV: Threshold Limit Value*

*PEL: Permissible Exposure Limit*

*REL: Recommended Exposure Limit*

*Ox. Liq. 2: Oxidizing liquids – Category 2*

*Met. Corr.1: Corrosive to metals – Category 1*

*Skin Corr. 1A: Skin corrosion/irritation – Category 1A*

*Skin Irrit. 2: Skin corrosion/irritation – Category 2*

*Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A*

· **Sources**

*Tables 3.1 and 3.2 from Annex 6 of EC 1272/2008, EC 1907/2006, EH40/2005 as amended 2011, Registry of Toxic Effects of Chemical Substances (RTECS), The Dictionary of Substances and their Effects, 1st Edition, IUCLID.*

- **Data compared to the previous version altered.** All sections have been updated.