

## **CERTIFICATE OF CONFORMANCE**

Agilent Product Name: BND Blue Scr Cap PTFE/Red Si Septa 100pk Agilent Part No: 5190-7024 Certificate Date: January 5, 2015

## **Quality Statement**

This product is certified and guaranteed to be fully compatible with all corresponding Agilent autosamplers and vials. All lots of this product have been manufactured under a Quality Management System registered to ISO 9001. (All corresponding Agilent glass vials are manufactured in a clean environment, using First Hydrolytic Type 1 borosilicate glass.) This product was manufactured to meet a number of critical specifications created by Agilent engineers to ensure a perfect fit with instrumentation.

## **Quality Test Results**

Critical Dimensions (Cap)	Cap OD	Pass
	Cap Width	Pass
	Septum Cavity ID	Pass
	Septum Retaining Ring ID	Pass
	Thru Hole ID	Pass
	Thread Major Diameter	Pass
	Thread Minor Diameter	Pass
Optoelectronic inspection	100%	Pass

Certified by:

heil househa

Neil Noronha Chemistries and Supplies Division Quality Manager