



CERTIFICATE OF CONFORMANCE

Agilent Product Name: Amb Rnd Bot HS Crimp Vial, 10mL – 100/pk

Agilent Part No: 5190-2238

Certificate Date: June 21, 2013

Quality Statement

This product is certified and guaranteed to be fully compatible with all corresponding Agilent autosamplers and closures. All lots of this product have been manufactured under a Quality Management System registered to ISO 9001. All glass vials are manufactured in a clean environment, using First Hydrolytic Type 1 borosilicate glass. This product was manufactured to meet a number of critical specifications created by Agilent engineers to ensure a perfect fit with instrumentation.

Quality Test Results

Critical Dimensions	Vial OD	Pass
	Vial Length	Pass
	Wall Weight	Pass
	Crimp OD (GPI – “T” Dim)	Pass
	Crimp Width (GPI – “A” Dim)	Pass
	Neck OD (GPI – “E” Dim)	Pass
	Neck ID (GPI – “I” Dim)	Pass
Optoelectronic inspection	100%	Pass

Certified by:

Neil Noronha
Chemistries and Supplies Division
Quality Manager