



## CERTIFICATE OF CONFORMANCE

**Agilent Product Name:** Screw caps for 6ml vials 100/PK

**Agilent Part No:** 9301-1379

**Certificate Date:** June 21, 2013

### Quality Statement

This product is certified and guaranteed to be fully compatible with all corresponding Agilent autosamplers and vials. All lots of this product have been manufactured under a Quality Management System registered to ISO 9001. (All corresponding Agilent glass vials are manufactured in a clean environment, using First Hydrolytic Type 1 borosilicate glass.) This product was manufactured to meet a number of critical specifications created by Agilent engineers to ensure a perfect fit with instrumentation.

### Quality Test Results

Critical Dimensions (Cap)	Cap OD	Pass
	Cap Width	Pass
	Septum Cavity ID	Pass
	Septum Retaining Ring ID	Pass
	Thru Hole ID	Pass
	Thread Major Diameter	Pass
	Thread Minor Diameter	Pass
Optoelectronic inspection	100%	Pass

**Certified by:**

Neil Noronha  
Chemistries and Supplies Division  
Quality Manager