



CERTIFICATE OF CONFORMANCE

Agilent Product Name: Cir Wide opening scrw top vial, 2ml 1000PK

Agilent Part No: 5183-2067

Certificate Date: June 21, 2013

Quality Statement

This product is certified and guaranteed to be fully compatible with all corresponding Agilent autosamplers and closures. All lots of this product have been manufactured under a Quality Management System registered to ISO 9001. All glass vials are manufactured in a clean environment, using First Hydrolytic Type 1 borosilicate glass. This product was manufactured to meet a number of critical specifications created by Agilent engineers to ensure a perfect fit with instrumentation.

Quality Test Results

| | | |
|---------------------------|-------------------------------------------------|------|
| Critical Dimensions | Vial OD | Pass |
| | Vial Length | Pass |
| | Interior Depth | Pass |
| | Wall Weight | Pass |
| | Thread Major Diameter (GPI – "T" Dim) | Pass |
| | Neck OD / Thread Minor Diameter (GPI – "E" Dim) | Pass |
| | Neck ID (GPI – "I" Dim) | Pass |
| Optoelectronic inspection | 100% | Pass |

Certified by:

Neil Noronha
Chemistries and Supplies Division
Quality Manager