

CHROMABOND® vacuum manifolds and accessories – compatibility chart

REF	Description	730150	730151	730360	730150N	730151N
730152	CHROMABOND® stainless steel needle, 12 pcs.	+	+	+	+	+
730153	CHROMABOND® collecting rack for manifold with 24 positions	-	+	-	-	+
730154	CHROMABOND® plastic needle PP, 12 pcs.	+	+	+	+	+
730157	CHROMABOND® collecting rack for manifold with 12 positions	+	-	-	-	-
730173	CHROMABOND® glass cabinet for manifold with 12 positions, w/o accessories	+	-	-	-	-
730174	CHROMABOND® glass cabinet for manifold with 16 or 24 positions, w/o accessories	-	-	+	-	-
730175	CHROMABOND® lid with gasket for manifold with 12 positions	+	-	-	-	-
730176	CHROMABOND® lid with gasket for manifold with 24 positions	-	+	-	-	+
730177	CHROMABOND® gaskets for lid for manifold with 12 positions	+	-	-	+	-
730178	CHROMABOND® gaskets for lid for manifold with 16 or 24 positions	-	+	+	-	+
730179	CHROMABOND® manometer with complete accessories	+	+	+	-	-
730183.12	CHROMABOND® Luer fittings female for lid, 12 pcs.	+	+	+	+	+
730184.12	CHROMABOND® Luer fittings male for lid, 12 pcs.	+	+	+	+	+
730185	CHROMABOND® valves plastic, 12 pcs.	+	+	+	+	+
730187	CHROMABOND® drying attachment with 12 positions	+	-	-	-	-
730188	CHROMABOND® drying attachment with 24 positions	-	+	-	-	+
730189.1	CHROMABOND® valve, brass, 1 pc.	+	+	+	+	+
730194	CHROMABOND® Luer caps for vacuum manifold, blue	+	+	+	+	+
730233	CHROMABOND® PP tank for 12 position vacuum manifold	+	-	-	-	-
730365	CHROMABOND® lid for manifold with 16 positions	-	-	+	-	-
730366	CHROMABOND® collecting rack for manifold with 16 positions	-	-	+	-	-
730564	CHROMABOND® PTFE connectors, 12 pcs.	+	+	+	+	+
730990	CHROMABOND® drying attachment for manifold with 16 positions	-	-	+	-	-
730106	CHROMABOND® stainless steel connector	+	+	+	+	+
730153N	CHROMABOND® collecting rack SPE for manifold with 24 positions	-	+	-	-	+
730154N	CHROMABOND® polypropylene needle SPE, 12 pcs.	+	+	+	+	+
730157N	CHROMABOND® collecting rack SPE for manifold with 12 positions	+	-	-	+	-
730173N	CHROMABOND® glass cabinet SPE for manifold with 12 positions, w/o accessories	-	-	-	+	-
730174N	CHROMABOND® glass cabinet SPE for manifold with 24 positions, w/o accessories	-	-	-	-	+
730175N	CHROMABOND® lid SPE with gasket for manifold with 12 positions	-	-	-	+	-
730176N	CHROMABOND® lid SPE with gasket for manifold with 24 positions	-	+	-	-	+
730177N	CHROMABOND® gaskets for lid SPE for manifold with 12 positions	+	-	-	+	-
730178N	CHROMABOND® gaskets for lid SPE for manifold with 24 positions	-	+	-	-	+
730179N	CHROMABOND® SPE manometer, complete	+	+	+	+	+
730185N	CHROMABOND® valves SPE, plastic, 12 pcs.	+	+	+	+	+
730187N	CHROMABOND® dry-cover SPE for manifold with 12 positions	-	-	-	+	-
730188N	CHROMABOND® dry-cover SPE for manifold with 24 positions	-	+	-	-	+
730233N	CHROMABOND® PP tank SPE for manifold with 12 positions, 2 pcs.	+	-	-	+	-
730564N	CHROMABOND® PP connectors SPE, 12 pcs.	+	+	+	+	+
730590N	CHROMABOND® PP tank SPE for manifold with 24 positions, 2 pcs.	-	+	-	-	+

+ compatible, - not compatible

730150 CHROMABOND® vacuum manifold for 12 columns
 730151 CHROMABOND® vacuum manifold for 24 columns
 730360 CHROMABOND® vacuum manifold for 16 columns
 730150N CHROMABOND® SPE vacuum manifold for 12 columns
 730151N CHROMABOND® SPE vacuum manifold for 24 columns

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG
 Valencienner Str. 11
 52355 Düren · Germany

DE Tel.: +49 24 21 969-0 info@mn-net.com
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com