


# NucleoSpin® Forensic Filters

*Mini spin solution for incubation and lysate separation of swabs, denim, cigarette butts,...*



**Certified forensic quality**  
EO treated plastic materials

**A convenient processing tool for solid sample carriers!**


- ▶ *Lysis and separation in one tube – no sample transfer, no extra pipetting steps*
- ▶ *Special filter design – absolutely leak-proof and compatible with all typical lysis buffers*
- ▶ *Complete DNA recovery and maximum process reliability – no DNA loss*
- ▶ *No cross-contamination due to collection tube with lid*
- ▶ *Ethylene oxide treated – free of amplifiable DNA*

**MACHERY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)



# NucleoSpin® Forensic Filters Procedure



Product (Mini spin columns)		Preps	REF
<b>NucleoSpin® Forensic Filters</b>	Mini spin filter for processing solid sample carriers	10 / 50 / 250	740988.10 / .50 / .250
Related products		Preps	REF
<b>NucleoSpin® Tissue</b>	Mini spin kit for isolation of genomic DNA from a broad range of samples	10 / 50 / 250	740952.10 / .50 / .250
<b>NucleoSpin® Tissue XS</b>	Mini spin kit for isolation of genomic DNA from extra small samples	10 / 50 / 250	740901.10 / .50 / .250
<b>NucleoSpin® 8 / 96 Trace</b>	8-well strips or 96-well plates for medium- / high-throughput DNA purification from forensic samples	12 x 8 / 60 x 8 2 x 96 / 4 x 96	740722.1 / .5 740726.2 / .4
<b>NucleoMag® Forensic</b>	Magnetic-bead based isolation of genomic DNA from forensic casework samples	1 x 96 / 4 x 96	744650.1 / .4
<b>NucleoMag® 96 Trace</b>	Magnetic-bead based isolation of genomic DNA from buccal swabs	1 x 96 / 4 x 96 / 24 x 96	744600.1 / .4 / .24

[www.mn-net.com](http://www.mn-net.com)

**MACHERY-NAGEL**



MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

**Germany**  
and international:  
Tel.: +49 24 21 969-0  
E-mail: [info@mn-net.com](mailto:info@mn-net.com)

**Switzerland:**  
MACHERY-NAGEL AG  
Tel.: +41 62 388 55 00  
E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

**France:**  
MACHERY-NAGEL EURL  
Tel.: +33 386 68 22 68  
E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

**USA:**  
MACHERY-NAGEL Inc.  
Tel.: +1 484 821 0984  
E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

