


Product Data Sheet

Date of Generation: 08.03.2016

Article-Nr.:	106 412
Type of Product:	HPLC desktop mount, electro conductive with gasket and locknut
Color:	black
Connection:	Ø32mm pipe, 1x 6,4- 9 mm ID (tube), 6x 2,3/3,2 OD (capillary)
Dimensions (W x D):	ca. 111 x Ø53 mm
Weight:	ca. 85 g
Material:	PE-HD
Surface Resistivity:	10 ⁴ - 10 ⁶ Ω



Post-/Lieferanschrift
S.C.A.T. Europe GmbH
Opelstraße 3
64546 Mörfelden
Deutschland

Telekommunikation
Tel.: +49-(0) 6105 - 305 586 - 0
Fax: +49-(0) 6105 - 305 586 - 99
e-Mail: info@scat-europe.com
Web: www.scat-europe.com

Bankverbindung
Volksbank Frankfurt
BLZ: 501 900 00
Konto-Nr.: 600 151 1678
SWIFT: FF VB DE FF

Gesellschaftsinformation
Sitz: 64546 Mörfelden
Amtsgericht: Darmstadt
Handelsregister B, Nr. 55045
Geschäftsführer: Johann Rittgasser