

# Datasheet

<b>Part No.:</b>	<b>106616</b>
<b>Description</b>	HPLC-Desktop mount with connecting thread GL25 (m), PE-HD electrostatic conductive, with 4x UNF 1/4" 28G and 2x NPT 1/8", pad diameter 53 mm, accessories must be ordered separately
<b>Measurements (WxDxH)</b>	54 x 54 x 90 mm
<b>Diameter:</b>	54,0 mm
<b>Weight:</b>	0,120 kg
<b>Colour</b>	black
<b>Material:</b>	PEHD-EL
<b>Chemical resistance:</b>	good
<b>Surface resistance:</b>	≤ 10 <sup>6</sup> Ohm
<b>Operating temperature:</b>	≤ 60° C
<b>Thread size:</b>	GL25

<b>Compliance / approvals</b>		
<b>Document No.:</b>	<b>Date</b>	<b>Inspection authority:</b>

<b>UN approval:</b>

<b>Other approvals / marks:</b>
<ul style="list-style-type: none"> <li>• seal with silicone</li> <li>• rotation lock: our recommendation for the boreholes, Ø 4-5mm diamond coring bit</li> </ul>

