

Injection port | Inlet liner PTV/LVI

Spot the difference in the SGE lineup!

Part number	Part description and detail
Agilent	
092210	Agilent AC Inlet Liners PTV Untreated, Length 73 mm PK5
092211	Agilent AC Inlet Liners PTV Untreated, Length 73 mm PK25
Agilent/Gerstel	
092243	Agilent Gerstel CIS 4 Inlet Liners LVI, Length 71 mm PK5
Bruker/Varian	
092245	Bruker/Varian 1078/1079 Inlet Liners LVI, Length 54 mm PK5
09224525	Bruker/Varian 1078/1079 Inlet Liners LVI, Length 54 mm PK25
PerkinElmer	
092097	PerkinElmer Inlet Liners PTV Gooseneck 1 mm ID, Length 88 mm PK5
PerkinElmer PSS	
092244	PerkinElmer PSS Inlet Liners PTV, Length 86.2 mm PK5
Shimadzu	
092042501	Shimadzu Inlet Liner PTV 3 Baffle, Length 95 mm PK1
0920425	Shimadzu Inlet Liners PTV 3 Baffle, Length 95 mm PK5
09229201	Shimadzu 2010 Inlet Liner LVI Straight Non-Sintered with Wool, Length 95 mm PK1
092292	Shimadzu 2010 Inlet Liners LVI Straight Non-Sintered with Wool, Length 95 mm PK5
Thermo Scientific TRACE 1300 Series GC, TRACE GC Ultra and FOCUS GC	
092155	Thermo Scientific Inlet Liners LVI, Length 120 mm PK5
09215525	Thermo Scientific Inlet Liners LVI, Length 120 mm PK25

Note - Bruker/Varian 1077 injectors use 78.5 mm "Agilent" liners

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com