

Separations | HPLC Column ProteCol[®] Amino

Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

Features and benefits

Allows basic compound separation in normal phase and carbohydrate analysis.

Recommended Applications

Reverse phase carbohydrate analysis and normal phase basic compound separations.

Product specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

Phase	Pore size (Å)	Particle size (µm)	Surface area	% Carbon load	USP code
Silica	120	5	350	4.0	L8

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

