

## Connections | GC ferrules graphite

For use on Ag	gilent and Thermo Scientific TRACE 1300 Series GC injectors and detectors at atmospheric pressure e.g. FID
072601	Graphite Short Ferrule with 1/4" Hole for 1/4" Fittings, Bruker/Varian 1075 & 1077 Injection Port Liners, Unions, & 1/4" OD Packed Columns, PK10
072602	Graphite Short Ferrule with 1/8" Hole for 1/8" Fittings, Unions, & 1/8" OD Packed Columns, PK10
072635	Graphite Short Reducing Ferrule with 0.5 mm Hole for Detector / Injectors (Not for MSD Interface Connections), PK10
072636	Graphite Short Reducing Ferrule with 0.8 mm Hole for Detector / Injectors (Not for MSD Interface Connections), PK10
For use on Br	ruker/Varian GC injectors and detectors at atmospheric pressure e.g. FID
072627	Graphite Ferrule with 0.5 mm Hole for 1/16" Fittings, PK10
072626	Graphite Ferrule with 0.8 mm Hole for 1/16" Fittings, PK10
)72622	Graphite Ferrule with 1/8" Hole for 1/8" OD Packed Columns, PK10
072621	Graphite Ferrule with 1/4" Hole for 1/4" OD Packed Columns, PK10
For use on St	nimadzu GC-2014, GC-2010, GC-17A and GC-14A injectors/detectors (not for MS interfaces or QP2010 injector)
0726080	Graphite Ferrule with 0.5 mm Hole for GC14A, GC17A, GC2010 and GC2014 Detector / Injectors
0726082	Graphite Ferrule with 0.8 mm Hole for GC14A, GC17A, GC2010 and GC2014 Detector / Injectors
0726001	Graphite Sealing Ring with 5 mm Hole for Shimadzu 14, 15A & 16 (Spl-14 Injector) Liners & 5 mm OD Packed Columns, PK10
or use on Th	nermo Scientific TRACE GC Ultra and FOCUS GC injectors and detectors at atmospheric pressure e.g. FID
072619	Graphite Ferrule with 0.5 mm Hole for Carlo Erba Detector / Injectors
072614	Graphite Ferrule with 0.8 mm Hole for Carlo Erba Detector / Injectors
)72622	Graphite Ferrule with 1/8" Hole for 1/8" OD Packed Columns, PK10
072621	Graphite Ferrule with 1/4" Hole for 1/4" OD Packed Columns, PK10
For use on Pe	erkinElmer injectors and detectors at atmospheric pressure e.g. FID
)72627	Graphite Ferrule with 0.5 mm Hole for 1/16" Fittings, PK10
)72624	Graphite Ferrules with 0.5 mm Hole For use with 1/8" Nuts on PerkinElmer Injectors & Detectors, PK10
072626	Graphite Ferrule with 0.8 mm Hole for 1/16" Fittings, PK10
0726280	Graphite Ferrules with 1.0 mm Hole For use with 1/8" Nuts on PerkinElmer Injectors & Detectors, PK10
072622	Graphite Ferrule with 1/8" Hole for 1/8" OD Packed Columns, PK10
072621	Graphite Ferrule with 1/4" Hole for 1/4" OD Packed Columns, PK10
Standard ferr	ule solutions - 1/16" fitting specification
072603	Graphite Short Ferrule with 1/16" Hole for 1/16" Fittings, Unions, PK10
072606	Graphite Short Reducing Ferrule with 1.2 mm Hole for 1/16" Fittings, PK10
072607	Graphite Short Reducing Ferrule with 1 mm Hole for 1/16" Fittings, PK10
072608	Graphite Short Reducing Ferrule with 0.8 mm Hole for 1/16" Fittings, PK10
072610	Graphite Short Reducing Ferrule with 0.5 mm Hole for 1/16" Fittings, PK10
072623	Graphite Ferrule with 1/16" Hole for 1/16" Fittings, PK10
072625	Graphite Ferrule with 1.0 mm Hole for 1/16" Fittings, PK10
Standard ferr	ule solutions - 1/8" fitting specification
0726270	Graphite Ferrule with 0.8 mm Hole for 1/8" Fittings, PK10
072628	Graphite Ferrule with 1.2 mm Hole for 1/8" Fittings, PK10
0726283	Graphite Ferrule with 1/16" Hole for 1/8" Fittings, PK10
072609	Graphite Short Reducing Ferrule with 1.2 mm Hole for 1/8" Fittings, PK10
	ule solutions - 1/4" fitting specification
072604	Graphite Short Reducing Ferrule with 6 mm Hole for 1/4" Fittings, PK10
0726151	Graphite Short Reducing Ferrule with 1/8" Hole for 1/4" Fittings, PK10
072617	Graphite Short Reducing Ferrule with 5/16" Hole for 1/4" Fittings, PK10
1726212	Graphite Reducing Ferrule with 6 mm Hole for 1/4" Fittings, PK 10
0726212 0726281	Graphite Reducing Ferrule with 6 mm Hole for 1/4" Fittings, PK10 Graphite Ferrule with 1/8" Hole for 1/4" Fittings, PK10

## Connections | GC ferrules graphite continued

Part number	Part description and detail
Other ferrules	
072629	Replacement Mini Union Graphite Ferrules with 0.8 mm Hole, PK10
072630	Replacement Mini Union Graphite Ferrules with 0.5 mm Hole, PK10
072605	Graphite Short Ferrule with 5 mm Hole for 5 mm Fittings, PK10
072613	Graphite Short Ferrule with 1/2" Hole for 1/2" Fittings, PK10
072615	Graphite Short Ferrule with 3/8" Hole for 3/8" Fittings, PK10
0726217	Graphite Sealing Ring for Varian 1078 & 1079 Injectors, PK10

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

