

Connections | HPLC Hexnut™

Quick and Easy HPLC Connections

Ideal for applications where corrosive solvents are being used, Hexnuts™ have inert contact surfaces, making them biocompatible. Stainless steel 10-32 thread fittings use a non-swaging Kel-F® or PEEK replaceable ferrule.



Features and benefits

- Stainless steel 10-32 thread fittings use a non-swaging Kei-F or PEEK replaceable ferrule.
- Inert contact surfaces make them bio compatible.
- Ideal for applications where corrosive solvents are being used.



Recommended applications

- Biotechnology.
- · Medical research.
- · Proteomics.
- Pharmaceuticals.
- Recommended for all HPLC connections, particularly where PEEK or PEEKsil® tubing is used.

Product specifications

- Maximum Pressure Rating: 5,000 psi
- Thread Type: 10-32

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com