

## Syringes | HPLC autosampler

Part number	Part description and detail	Replacement needle (*needle and plunger kit)	Replacement plunger
<b>CTC Analytics</b>			
002710	10 µL fixed needle CTC syringe with 5.1 cm 0.72 mm OD LC needle	-	-
003715	25 µL fixed needle CTC syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle (8 mm OD barrel)	-	0318922
004810	50 µL fixed needle CTC syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	-	031821
005715	100 µL fixed needle CTC syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	-	0318261
005720	100 µL fixed needle CTC syringe with GT plunger and 5.1 cm 0.72 mm OD (0.4 mm ID) LC needle	-	0318261
006330	250 µL fixed needle CTC syringe with GT plunger and 5.1 cm 0.72 mm OD (0.4 mm ID) LC needle	038250	031831
007720	500 µL fixed needle CTC syringe with GT plunger and 5.1 cm 0.72 mm OD (0.4 mm ID) LC needle	-	0318928
<b>CTC Analytics and Thermo Scientific</b>			
005331	100 µL fixed needle CTC/Thermo syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	-	0318261
005330	100 µL removable needle CTC/Thermo syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	038250	031826
<b>CTC Analytics PAL RTC and Thermo Scientific TriPlus RSH</b>			
002880	10 µL removable needle CTC RTC and Thermo RSH syringe with GT plunger and 5.7 cm 0.72 mm OD LC needle	037880	0318121
008150	1 mL fixed needle CTC RTC and Thermo RSH syringe with GT plunger and 5.7 cm 0.72 mm OD LC tipped needle	-	032842
008151	1 mL removable needle CTC RTC and Thermo RSH syringe with GT plunger and 5.7 cm 0.72 mm OD LC tipped needle	039280	032842
008950	10 mL fixed needle CTC RTC and Thermo RSH syringe with GT plunger and 5.7 cm 1.07 mm OD LC tipped needle	-	032855
008951	10 mL removable needle CTC RTC and Thermo RSH syringe with GT plunger and 5.7 cm 1.07 mm OD LC tipped needle	031590	032855
<b>Hitachi</b>			
007660	500 µL syringe with GT plunger and M10 termination	-	031837
<b>PerkinElmer</b>			
004995	50 µL syringe with GT plunger and 1/4-28 UNF termination	-	-
005990	100 µL syringe with GT plunger and 1/4-28 UNF termination	-	-
006995	250 µL syringe with GT plunger and 1/4-28 UNF termination	-	031833
007995	500 µL syringe with GT plunger and 1/4-28 UNF termination	-	-
008185	1 mL syringe with GT plunger and 1/4-28 UNF termination	-	0318441
<b>Spark</b>			
006683	250 µL Spark syringe with GT plunger	-	0318348
<b>Thermospec</b>			
006660	250 µL Thermospec syringe	-	-
<b>Waters</b>			
003990	25 µL Waters WISP syringe with GT plunger	-	031819
006690	250 µL Waters WISP syringe with GT plunger	-	031834

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

*Specifications are subject to change without notice.*

*Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.*