

## Syringes | High pressure

Part number	Part description and detail	Replacement needle
008170	1 mL removable needle high pressure gas tight syringe with gas tight ON/OFF valve and 7 cm 0.63 mm OD bevel tipped needle	0315720
008171	1 mL removable needle high pressure gas tight syringe with gas tight ON/OFF valve and 7 cm 0.63 mm OD sleeved side hole needle	0315722

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

*Specifications are subject to change without notice.*

*Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.*