

Injection port | Inlet liner straight

Spot the difference in the SGE lineup!

| Part number | Part description and detail |
|--|---|
| Agilent, Bruker/Varian 1077 and Thermo Scientific TRACE 1300 series | |
| 09201601 | Inlet Liner 1.2 mm ID Straight, Length 78.5 mm PK1 |
| 092016 | Inlet Liners 1.2 mm ID Straight, Length 78.5 mm PK5 |
| 092224 | Inlet Liners 1.2 mm ID Straight, Length 78.5 mm PK25 |
| 09200401 | Inlet Liner, Quartz, 2 mm ID Straight, Length 78.5 mm PK1 |
| 092004 | Inlet Liners, Quartz, 2 mm ID Straight, Length 78.5 mm PK5 |
| 09200701 | Inlet Liner 4 mm ID Straight, Length 78.5 mm PK1 |
| 092007 | Inlet Liners 4 mm ID Straight, Length 78.5 mm PK5 |
| 092222 | Inlet Liners 4 mm ID Straight, Length 78.5 mm PK25 |
| 09200101 | Inlet Liner with Quartz Wool 4 mm ID, Length 78.5 mm PK1 |
| 092001 | Inlet Liners with Quartz Wool 4 mm ID, Length 78.5 mm PK5 |
| 092220 | Inlet Liners with Quartz Wool 4 mm ID, Length 78.5 mm PK25 |
| Agilent | |
| 09223401 | Agilent Inlet Liner Packed Column, Length 92 mm PK1 |
| 092234 | Agilent Inlet Liners Packed Column, Length 92 mm PK5 |
| 092235 | Agilent Inlet Liners Packed Column, Length 92 mm PK25 |
| Bruker/Varian | |
| 092031 | Bruker/Varian 1078/1079 Inlet Liners 0.5 mm ID Straight, Length 54 mm PK5 |
| 092117 | Bruker/Varian 1078/1079 Inlet Liners SPME 0.75 mm ID, Length 54 mm PK5 |
| 092027 | Bruker/Varian 1093/1094 Inlet Liners 0.5 mm ID, Length 54 mm PK5 |
| 092030 | Bruker/Varian 1093/1094 Inlet Liners 0.8 mm ID (rest 0.25mm), Length 54 mm PK5 |
| 092034 | Bruker/Varian 1093/1094 Inlet Liners 0.8 mm ID (rest 0.50mm), Length 54 mm PK5 |
| Dani | |
| 092500 | Dani Inlet Liners with Quartz Wool 4 mm ID, Length 70 mm PK5 |
| 092501 | Dani Inlet Liners 4 mm ID Straight, Length 70 mm PK5 |
| PerkinElmer | |
| 092098 | PerkinElmer PSS Inlet Liners 2 mm ID, Length 86.2 mm PK5 |
| 09210301 | PerkinElmer Inlet Liner 2 mm ID Straight, Length 92 mm PK1 |
| 092103 | PerkinElmer Inlet Liners 2 mm ID Straight, Length 92 mm PK5 |
| 09210001 | PerkinElmer Inlet Liner 4 mm ID Straight, Length 92 mm PK1 |
| 092100 | PerkinElmer Inlet Liners 4 mm ID Straight, Length 92 mm PK5 |
| 09210025 | PerkinElmer Inlet Liners 4 mm ID Straight, Length 92 mm PK25 |
| 092236 | PerkinElmer Inlet Liners Packed Column, Length 112 mm PK5 |
| Shimadzu | |
| 09208401 | Shimadzu 9A/16A Inlet Liner 3.4 mm ID, Length 139 mm PK1 |
| 092084 | Shimadzu 9A/16A Inlet Liners 3.4 mm ID, Length 139 mm PK5 |
| 09208901 | Shimadzu 2010/17A Inlet Liner 0.75 mm ID, Length 95 mm PK1 |
| 092089 | Shimadzu 2010/17A Inlet Liners 0.75 mm ID, Length 95 mm PK5 |
| 092086101 | Shimadzu 2010/17A Inlet Liner 2.6 mm ID, Length 95 mm PK1 |
| 0920861 | Shimadzu 2010/17A Inlet Liners 2.6 mm ID, Length 95 mm PK5 |
| 09206401 | Shimadzu 2010/17A Inlet Liner 3.4 mm ID, Length 95 mm PK1 |
| 092064 | Shimadzu 2010/17A Inlet Liners 3.4 mm ID, Length 95 mm PK5 |
| 09229201 | Shimadzu 2010/17A Inlet Liner Non-Sintered LVI with Wool Straight, Length 95 mm PK1 |

Injection port | Inlet liner straight continued

| Part number | Part description and detail |
|--|--|
| Thermo Scientific TRACE GC Ultra and FOCUS GC | |
| 09214801 | Thermo Scientific Inlet Liner 0.8 mm ID Straight, Length 105 mm PK1 |
| 092148 | Thermo Scientific Inlet Liners 0.8 mm ID Straight, Length 105 mm PK5 |
| 092147 | Thermo Scientific Inlet Liners 3 mm ID, Length 105 mm PK5 |
| 09215001 | Thermo Scientific Inlet Liner 5 mm ID Straight, Length 105 mm PK1 |
| 092150 | Thermo Scientific Inlet Liners 5 mm ID Straight, Length 105 mm PK5 |
| 092142 | Thermo Scientific TRACE Inlet Liners 2 mm ID, Length 120 mm PK5 |

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com