

Injection port | Inlet liner tapered/gooseneck

Spot the difference in the SGE lineup!

Part number	Part description and detail
Agilent and Thermo Scientific TRACE 1300 GC	
092013	Inlet Liners Gooseneck 2 mm ID, Length 78.5 mm PK5
09201001	Inlet Liner Gooseneck with Quartz Wool 4 mm ID, Length 78.5 mm PK1
092010	Inlet Liners Gooseneck with Quartz Wool 4 mm ID, Length 78.5 mm PK5
092223	Inlet Liners Gooseneck with Quartz Wool 4 mm ID, Length 78.5 mm PK25
09201701	Inlet Liner Tapered 4 mm ID, Length 78.5 mm PK1
092017	Inlet Liners Tapered 4 mm ID, Length 78.5 mm PK5
092229	Inlet Liners Tapered 4 mm ID, Length 78.5 mm PK25
09201901	Inlet Liner Tapered with Quartz Wool 4 mm ID, Length 78.5 mm PK1
092019	Inlet Liners Tapered with Quartz Wool 4 mm ID, Length 78.5 mm PK5
092218	Inlet Liners Tapered with Quartz Wool 4 mm ID, Length 78.5 mm PK25
0926201901	Inlet Liner Tapered with Quartz Wool Pre-fitted O-ring and TI Deactivation 4 mm ID, Length 78.5 mm PK1
09262019	Inlet Liners Tapered with Quartz Wool Pre-fitted O-ring and TI Deactivation 4 mm ID, Length 78.5 mm PK5
Bruker/Varian	
092024	Bruker/Varian 1075/1077 Inlet Liners Restricted 2 mm ID, Length 74 mm PK5
092228	Bruker/Varian 1075/1077 Inlet Liners Restricted 2 mm ID, Length 74 mm PK25
09202101	Bruker/Varian 1075/1077 Inlet Liner Restricted with Quartz Wool 4 mm ID, Length 72 mm PK1
092021	Bruker/Varian 1075/1077 Inlet Liners Restricted with Quartz Wool 4 mm ID, Length 72 mm PK5
09203801	Bruker/Varian 1078/1079 Inlet Liner Tapered 3.4 mm ID, Length 54 mm PK1
092038	Bruker/Varian 1078/1079 Inlet Liners Tapered 3.4 mm ID, Length 54 mm PK5
09203825	Bruker/Varian 1078/1079 Inlet Liners Tapered 3.4 mm ID, Length 54 mm PK25
09203901	Bruker/Varian 1078/1079 Inlet Liner Tapered 2 mm ID, Length 54 mm PK1
092039	Bruker/Varian 1078/1079 Inlet Liners Tapered 2 mm ID, Length 54 mm PK5
Dani	
092504	Dani Inlet Liners Tapered 4 mm ID, Length 70 mm PK5
Optic	
09227301	Optic Inlet Liner with Frit 3 mm ID, Length 81 mm PK1
092273	Optic Inlet Liners with Frit 3 mm ID, Length 81 mm PK5
PerkinElmer	
09209701	PerkinElmer Inlet Liner PTV Gooseneck 1 mm ID, Length 88 mm PK1
092097	PerkinElmer Inlet Liners PTV Gooseneck 1 mm ID, Length 88 mm PK5
092091	PerkinElmer Inlet Liners Gooseneck 3 mm ID, Length 100 mm PK5 (for 5000 and Sigma Series)
0920990	PerkinElmer Inlet Liners Tapered 4 mm ID, Length 92 mm PK5
Shimadzu	
09208701	Shimadzu 2010/17A Inlet Liner Tapered 2.6 mm ID, Length 95 mm PK1
092087	Shimadzu 2010/17A Inlet Liners Tapered 2.6 mm ID, Length 95 mm PK5
09208501	Shimadzu 2010/17A Inlet Liner Gooseneck 3.4 mm ID, Length 95 mm PK1
092085	Shimadzu 2010/17A Inlet Liners Gooseneck 3.4 mm ID, Length 95 mm PK5
09208525	Shimadzu 2010/17A Inlet Liners Gooseneck 3.4 mm ID, Length 95 mm PK25
09206101	Shimadzu 2010/17A Inlet Liner Gooseneck with Quartz Wool 3.4 mm ID, Length 95 mm PK1
092061	Shimadzu 2010/17A Inlet Liners Gooseneck with Quartz Wool 3.4 mm ID, Length 95 mm PK5
09207101	Shimadzu 2010/17A Inlet Liner Tapered 3.4 mm ID, Length 95 mm PK1
092071	Shimadzu 2010/17A Inlet Liners Tapered 3.4 mm ID, Length 95 mm PK5
092077	Shimadzu 2010/17A Inlet Liners Gooseneck 3.4 mm ID, Length 95 mm PK5
0920831	Shimadzu 14 Inlet Liners Tapered 3.4 mm ID, Length 99 mm PK5
092083125	Shimadzu 14 Inlet Liners Tapered 3.4 mm ID, Length 99 mm PK25
092082	Shimadzu 14 Inlet Liners Mid-Gooseneck 3.4 mm ID, Length 99 mm PK5

Injection port | Inlet liner tapered/gooseneck continued

Part number	Part description and detail
Thermo Scientific TRACE GC Ultra and FOCUS GC	
09214125	Thermo Scientific Inlet Liners Tapered 3 mm ID, Length 105 mm PK25
09214401	Thermo Scientific Inlet Liner Tapered 5 mm ID, Length 105 mm PK1
092144	Thermo Scientific Inlet Liners Tapered 5 mm ID, Length 105 mm PK5
09214425	Thermo Scientific Inlet Liners Tapered 5 mm ID, Length 105 mm PK25

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com