

## Syringes | HPLC autosampler

Known for 50 years in chromatography syringes, the SGE range of HPLC autosampler syringes are designed and tested to meet critical autosampler specifications.

## **Product specifications**

Available for CTC Analytics, Spark Holland, Hitachi, Thermo Scientific, PerkinElmer, Waters WISP, HPLC autosamplers.

Accuracy and reproducibility	±1% (dispensed volume)
Temperature range	5-70°C fixed needle 5-120°C removable needle

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.

Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.



