

# Syringes | MS

#### The choice for sensitivity

Understanding how important your MS results are, and how precious your sample is, the SGE MS syringe and needle range is designed to improve the sensitivity of your MS results.

#### Features and benefits

- Lowest syringe carryover.
- · Reduced non specific sample interaction and increased analyte recovery.
- Ultra smooth removable needles with choice of surface treatment.

## Recommended applications

MS syringes and needles are MS ready.

### **Product specifications**

Available in 10 µL to 500 µL manual syringes, 10 µL to 100 µL CTC Analytics and Thermo Scientific autosampler syringes, 10 µL and 25 µL Agilent autosampler syringes, 10 µL Shimadzu autosampler syringes, ultra smooth, ultra smooth hydrophilic, and ultra smooth hydrophobic needles.

Accuracy and reproducibility	±1% (dispensed volume)
Temperature range	5-120°C
International standards traceability	

# For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.

Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.



