

Tubing | MyCapLC™ Kit

Customize CapLC for Your Application

MyCapLC™ Kit enables you to create your own capillary columns with your HPLC pump and the simple packing instructions. You can easily pack your own capillary columns in a matter of hours.

Features and benefits

- DIY capillary column packing, no special pumps or instrumentation required.
- Quick, convenient and cost effective solution for capillary columns.
- Robust and easy to pack hardware.

Recommended Applications

Perfect for life sciences and omic analysis.

Product specifications

MyCapLC Kit (suitable to pack three columns):

- 0.1 g C18 3 μm 200 Å
- 1 x Slurry reservoir (50 mm x 2.1 mm ID column tubing)
- 2 x End fittings with inserted PEEK sealing rings
- 1 x PEEK Fingertight fitting
- 3 x Column bodies (1/16" OD PEEKsil®, lengths from 50 mm- 150 mm and IDs from 150 μm - 530 μm)
- 6 x Connection capillaries (1/32" and 1/16" OD options available in 50 μm ID Inlet length is 250 mm Outlet length is 100 mm)
- 6 x Male unions (1/16")
- 6 x Female unions (1/32")
- 6 x PEEK double ended ferrules (1/16"-1/32") with 0.5 μm woven mesh filters
- 2 x Wrenches

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

