

## Accessories | Needle cleaning kits

The SGE syringe is a precision instrument with a high standard of fit between the plunger and the glass barrel. Like most precision instruments, regular maintenance is important for ensuring long life and robust performance.

The SGE syringe needle cleaning kit contains a range of stylet wires for needle cleaning, tweezers and a surfactant for barrel cleaning.



### Features and benefits

Convenient package with everything needed to clean SGE syringe needles.

### Recommended applications

The needle cleaning kit is ideal for cleaning and unblocking needles.

### Product specifications

A range of stylet wires, tweezers for holding the wires and a phosphate-free, non-ionic surfactant are provided for cleaning syringe needles.

#### Kit includes:

- 15 cm long x 0.10 mm OD wires, PK5
- 15 cm long x 0.15 mm OD wires, PK5
- 15 cm long x 0.20 mm OD wires, PK5
- 20 mL cleaning solution concentrate
- Tweezers

Replacement stylet wires are available.

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

*Specifications are subject to change without notice.*

*Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.*