

Separations | HPLC Guard Column ProteCol®

Turn imagination into results

Part number	Part description and detail
Stainless Steel	
C8H	
2C83-05G02S	C8H x 5 µm Particle Size x 120 Å Pore Size, 10 mm Length x 2.1 mm ID, Stainless Steel ProteCol HPLC Guard Column PK3
C18G	
2C185-03G04S	C18G x 3 µm Particle Size x 120 Å Pore Size, 10 mm Length x 4.0 mm ID, Stainless Steel ProteCol HPLC Guard Column PK3
2C185-05G04S	C18G x 5 µm Particle Size x 120 Å Pore Size, 10 mm Length x 4.0 mm ID, Stainless Steel ProteCol HPLC Guard Column PK3
2C181-03G02S	C18G x 3 µm Particle Size x 200 Å Pore Size, 10 mm Length x 2.1 mm ID, Stainless Steel ProteCol HPLC Guard Column PK3
C18H	
2C182-03G02S	C18H x 3 µm Particle Size x 120 Å Pore Size, 10 mm Length x 2.1 mm ID, Stainless Steel ProteCol HPLC Guard Column PK3
2C182-05G04S	C18H x 5 µm Particle Size x 120 Å Pore Size, 10 mm Length x 4.0 mm ID, Stainless Steel ProteCol HPLC Guard Column PK3
C18Q	
2C183-05G04S	C18Q x 5 µm Particle Size x 100 Å Pore Size, 10 mm Length x 4.0 mm ID, Stainless Steel ProteCol HPLC Guard Column PK3
Silica	
2SIL-05G04S	Silica x 5 µm Particle Size x 120 Å Pore Size, 10 mm Length x 4.0 mm ID, Stainless Steel ProteCol HPLC Guard Column PK3
PEEKsil	
C18G	
2C184-03G30K	C18G x 3 µm Particle Size x 300 Å Pore Size, 10 mm Length x 300 µm ID, PEEKsil ProteCol HPLC Guard Column PK3

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com