

Connections | SilTite ferrules

Part number	Part description and detail
For use on Varian/Bruker interface connections	
073300	SilTite metal - initial installation 0.10-0.25mm ID columns for Varian/Bruker GCMS interface connections (10 ferrules, 2 nuts)
073301	SilTite metal - initial installation 0.32 mm ID columns for Varian/Bruker GCMS interface connections (10 ferrules, 2 nuts)
073302	SilTite metal - initial installation 0.45-0.53 mm ID columns for Varian/Bruker GCMS interface connections (10 ferrules, 2 nuts)
073231	SilTite metal nuts for Varian/Bruker GCMS interface connections, PK5
For use on Shimadzu	
073204	SilTite metal - initial installation for Shimadzu QP5000/5050 standard MS interface for use with 0.10-0.25 mm ID columns (10 ferrules, 2 nuts)
073205	SilTite metal - initial installation for Shimadzu QP5000/5050 standard MS interface for use with 0.32 mm ID columns (10 ferrules, 2 nuts)
073233	SilTite metal nuts - for use with QP5000/5050 standard MS interface, PK5
073227	SilTite ferrules for Shimadzu GC-2014, GC-2010, GC-17A and GC-14A detector/injectors, and QP5000/5050 standard MS interface (not for MS interfaces or QP2010 injector) for use with 0.10-0.25 mm ID columns, PK10
073228	SilTite ferrules for Shimadzu GC-2014, GC-2010, GC-17A and GC-14A detector/injectors, and QP5000/5050 standard MS interface (not for MS interfaces or QP2010 injector) for use with 0.32 mm ID columns, PK10
073229	SilTite ferrules for Shimadzu GC-2014, GC-2010, GC-17A and GC-14A detector/injectors, and QP5000/5050 standard MS interface (not for MS interfaces or QP2010 injector) for use with 0.45-0.53 mm ID columns, PK10
073350	SilTite metal - initial installation for Shimadzu GC-2014, GC-2010, GC-17A and GC-14A detector/injectors, and QP5000/5050 standard MS interface (not for MS interfaces or QP2010 injector) for 0.10-0.25 mm ID columns (10 ferrules, 2 nuts)
073351	SilTite metal - initial installation for Shimadzu GC-2014, GC-2010, GC-17A and GC-14A detector/injectors, and QP5000/5050 standard MS interface (not for MS interfaces or QP2010 injector) for 0.32 mm ID columns (10 ferrules, 2 nuts)
073352	SilTite metal - initial installation for Shimadzu GC-2014, GC-2010, GC-17A and GC-14A detector/injectors, and QP5000/5050 standard MS interface (not for MS interfaces or QP2010 injector) for 0.45-0.53 mm ID columns (10 ferrules, 2 nuts)
073232	SilTite metal nuts, slotted, for Shimadzu GC-2014, GC-2010, GC-17A and GC-14A detector/injectors, and QP5000/5050 standard MS interface (not for MS interfaces or QP2010 injector), PK5
For use on Thermo Scientific interface connections	
073200	SilTite metal - initial installation for GCMS interface connection for use with 0.10-0.25 mm ID columns (10 ferrules, 2 nuts)
073201	SilTite metal - initial installation for GCMS interface connection for use with 0.32 mm ID columns (10 ferrules, 2 nuts)
073202	SilTite metal - initial installation for GCMS interface connection for use with 0.45-0.53 mm ID columns, PK10 (10 ferrules, 2 nuts)
073270	SilTite metal - initial installation for Thermo Scientific TRACE 1300 series GC split/splitless injectors for use with 0.10-0.25 mm ID columns (10 ferrules, 2 nuts, 2 inlet base seals)
073271	SilTite metal - initial installation for Thermo Scientific TRACE 1300 series GC split/splitless injectors for use with 0.32 mm ID columns (10 ferrules, 2 nuts, 2 inlet base seals)
073272	SilTite metal - initial installation for Thermo Scientific TRACE 1300 series GC split/splitless injectors for use with 0.45-0.53 mm ID columns (10 ferrules, 2 nuts, 2 inlet base seals)
073273	SilTite metal - initial installation for Thermo Scientific TRACE 1300 series GC split/splitless injectors for use with 1/32" columns (10 ferrules, 2 nuts, 2 inlet base seals)
073220	Replacement SilTite metal ferrules for all Thermo Scientific GCMS interface connections for use with 0.10-0.25 mm ID columns, PK10
073221	Replacement SilTite metal ferrules for all Thermo Scientific GCMS interface connections for use with 0.32 mm ID columns, PK10
073222	Replacement SilTite metal ferrules for all Thermo Scientific GCMS interface connections for use with 0.45-0.53 mm ID columns, PK10
073219	Replacement SilTite metal ferrules for all Thermo Scientific GCMS interface connections for use with 1/32" columns, PK10
073220	Replacement SilTite metal ferrules for Thermo Scientific TRACE 1300 series GC split/splitless injectors for use with 0.10-0.25 mm ID columns, PK10
073221	Replacement SilTite metal ferrules for Thermo Scientific TRACE 1300 series GC split/splitless injectors for use with 0.32 mm ID columns, PK10
073222	Replacement SilTite metal ferrules for Thermo Scientific TRACE 1300 series GC split/splitless injectors for use with 0.45-0.53 mm ID columns, PK10
073219	Replacement SilTite metal ferrules for Thermo Scientific TRACE 1300 series GC split/splitless injectors for use with 1/32" columns, PK10

Part number	Part description and detail
073224	Replacement SilTite metal nuts for all GCMS interface connections, PK5
073226	Replacement SilTite metal nuts for Thermo Scientific TRACE 1300 series GC split/splitless injector, PK5
073400	Replacement SilTite base seals For Thermo Scientific TRACE 1300 series GC split/splitless injector, PK2
073401	Replacement SilTite base seals For Thermo Scientific TRACE 1300 series GC split/splitless injector, PK10
For use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface	
073219	SilTite ferrules for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface for use with 1/32" columns, PK10
073220	SilTite ferrules for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface for use with 0.10-0.25 mm ID columns, PK10
073221	SilTite ferrules for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface for use with 0.32 mm ID columns, PK10
073222	SilTite ferrules for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface for use with 0.45-0.53 mm ID columns, PK10
073223	SilTite ferrules 1/16", for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface, PK10
073200	SilTite metal - initial installation for 0.10-0.25 mm ID columns (10 ferrules, 2 nuts)
073201	SilTite metal - initial installation for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface for use with 0.32 mm ID columns (10 ferrules, 2 nuts)
073202	SilTite metal - initial installation for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface for use with 0.45-0.53 mm ID columns (10 ferrules, 2 nuts)
073203	SilTite metal - initial installation 1/16" for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface (10 ferrules, 2 nuts)
073224	SilTite metal nuts 10/32" 0.8 mm ID for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface, PK5
073225	SilTite metal nuts 10/32" 1/16" ID for use on PerkinElmer MS interface, Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface, Agilent MS interface, PK5

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com