

Accessories | Syringe valves

Extraction to injection in a single process

Part number	Part description and detail
Syringe valves	
031905	Push-pull valve with replaceable 5 cm 0.63 mm OD cone tip needle for 25 μ L to 2.5 mL removable needle syringes
031907	Push-pull valve with Luer Lock termination for 25 μ L to 2.5 mL removable needle syringes
031910	Push-button valve with Luer Lock termination for 5 mL to 2 L removable Luer Lock or removable needle syringes
031915	Push-button valve with Luer Lock termination for any Luer Lock and Luer tip syringes or other Luer devices
Replacement needles	
038810	5 cm 0.63 mm OD cone tipped needle for valve 031905
038815	5 cm 0.63 mm OD side hole dome tipped needle for valve 031905

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.

Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.

