

## Syringes | Headspace soil gas tight

Part number	Part description and detail	Replacement needle	Replacement plunger	Replacement probe	Replacement valve
008775	5 mL headspace soil gas tight syringe with GT plunger and push-button valve with Luer Lock termination, removable side hole Luer probe, and 5 cm 0.63 mm OD bevel tipped needle	039802	031856	031571	031910
008975	10 mL headspace soil gas tight syringe with GT plunger and push-button valve with Luer Lock termination, removable side hole Luer probe, and 5 cm 0.63 mm OD bevel tipped needle	039802	031862	031571	031910
<b>Replacement parts</b>					
039802	5 cm 0.63 mm OD bevel tipped needle for Luer Lock, PK5				
031856	Replacement GT plunger for 5 mL headspace soil gas tight syringe 008775				
031862	Replacement GT plunger for 10 mL headspace soil gas tight syringe 008975				
031571	6 cm side hole Luer probe for headspace soil gas tight syringe				
031910	Push-button valve with Luer Lock termination for 5 mL to 2 L removable Luer Lock or removable needle syringes				

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

*Specifications are subject to change without notice.*

*Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.*