

Syringes | For Rheodyne® and Valco® Valves

Part number	Part description and detail	Replacement needle (*needle and plunger kit)	Replacement plunger
001301	5 µL fixed needle syringe with 5.1 cm 0.72 mm OD LC needle	-	-
002301	10 µL fixed needle syringe with 5.1 cm 0.72 mm OD LC needle	-	-
002315	10 µL fixed needle syringe with 5.1 cm 0.72 mm OD LC needle PK6	-	-
002300	10 µL fixed needle syringe with SuperFlex plunger and 5.1 cm 0.72 mm OD LC needle	-	-
002330	10 µL fixed needle syringe with SuperFlex plunger and 5.1 cm 0.72 mm OD LC needle PK6	-	-
002335	10 µL fixed needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	-	031801
002313	10 µL removable needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	037250	031811
003300	25 µL fixed needle syringe with 5.1 cm 0.72 mm OD LC needle	-	-
003312	25 µL removable needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	038250	031815
004300	50 µL fixed needle syringe with 5.1 cm 0.72 mm OD LC needle	-	-
004312	50 µL removable needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	038250	031820
005300	100 µL fixed needle syringe with 5.1 cm 0.72 mm OD LC needle	-	-
005312	100 µL removable needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	038250	031825
006300	250 µL fixed needle syringe with 5.1 cm 0.72 mm OD LC needle	-	-
006310	250 µL removable needle syringe with 5.1 cm 0.72 mm OD LC needle	038250	-
006312	250 µL removable needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	038250	031830
007300	500 µL fixed needle syringe with and 5.1 cm 0.72 mm OD LC needle	-	-
007310	500 µL removable needle syringe with 5.1 cm 0.72 mm OD LC needle	038250	-
007312	500 µL removable needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	038250	031835
008105	1 mL removable needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	039250	031842
008505	2.5 mL removable needle syringe with GT plunger and 5.1 cm 0.72 mm OD LC needle	039250	031852

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.

Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.