

# Syringes | Syringes with valves

Extraction to injection in a single process

Part number	Part description and detail	Replacement needle	Replacement plunger	Replacement valve
<b>Syringes pre-fitted with push-pull valve with replaceable needle</b>				
004279	50 µL syringe with GT plunger and push-pull valve with 5 cm 0.63 mm OD cone tipped needle	038810	031820	031905
005279	100 µL syringe with GT plunger and push-pull valve with 5 cm 0.63 mm OD cone tipped needle	038810	031825	031905
006279	250 µL syringe with GT plunger and push-pull valve with 5 cm 0.63 mm OD cone tipped needle	038810	031830	031905
007279	500 µL syringe with GT plunger and push-pull valve with 5 cm 0.63 mm OD cone tipped needle	038810	031835	031905
008110	1 mL syringe with GT plunger and push-pull valve with 5 cm 0.63 mm OD cone tipped needle	038810	031842	031905
008510	2.5 mL syringe with GT plunger and push-pull valve with 5 cm 0.63 mm OD cone tipped needle	038810	031852	031905
<b>Syringes pre-fitted with Luer Lock valve</b>				
008160	1 mL syringe with GT plunger push-pull valve with Luer Lock termination		031842	031907
008560	2.5 mL syringe with GT plunger push-pull valve with Luer Lock termination		031852	031907
008770	5 mL syringe with GT plunger and push-button valve with Luer Lock termination		031856	031910
008970	10 mL syringe with GT plunger and push-button valve with Luer Lock termination		031862	031910
009472	25 mL syringe with GT plunger and push-button valve with Luer Lock termination		031870	031910
009670	50 mL syringe with GT plunger and push-button valve with Luer Lock termination		0312170	031910
009770	100 mL syringe with GT plunger and push-button valve with Luer Lock termination		0312176	031910
<b>Replacement valves</b>				
031905	Push-pull valve with replaceable 5 cm 0.63 mm OD cone tip needle for 25 µL to 2.5 mL removable needle syringes			
031907	Push-pull valve with Luer Lock termination for 25 µL to 2.5 mL removable needle syringes			
031910	Push-button valve with Luer Lock termination for 5 mL to 2 L removable Luer Lock or removable needle syringes			
031915	Push-button valve with Luer Lock termination for any Luer Lock and Luer tip syringes or other Luer devices			
<b>Replacement needles</b>				
038810	5 cm 0.63 mm OD cone tipped needle for valve 031905			
038815	5 cm 0.63 mm OD side hole dome tipped needle for valve 031905			

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

*Specifications are subject to change without notice.*

*Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.*