

Ion Transfer Tube

Fusionシリーズ / TSQ Quantiva / TSQ Endula

Instrument	P/N	ID (mm)	Length (mm)	Seal/o-ring
Orbitrap Fusion Lumos TSQ Quantiva	79005-20641	0.58 x 1.25	58.42	70005-20922
Orbitrap Fusion TSQ Endura	70005-20606	0.58	58.42	97055-20442

LTQ シリーズ / Orbitrap シリーズ / Exactiveシリーズ

Instrument	P/N	ID (mm)	Length (mm)
Orbitrap Elite / Orbitrap Velos Pro / Orbitrap Velos / Velos Pro / Velos / Q Exactive HF / Q Exactive Plus / Q Exactive / Q Exactive Focus / Exactive Plus	70005-20606	0.58	58.42
LTQ Orbitrap XL / LTQ Orbitrap Discovery / LTQ XL	97055-20198	0.55	101.6
LCQ Fleet / LXQ	97055-20517	0.50	123.7

TSQ Quantumシリーズ

Instrument	P/N	ID (mm)	Length (mm)
TSQ Vantage	70005-20423	0.58	58.42
TSQ Quantum Access Max	70111-20972	0.50	101.6
TSQ Quantum Ultra	97055-20199	0.58	101.6
TSQ Quantum Discovery Max	70111-20396	0.50	82.55

